

covetrus ™

Covetrus Branded Products

2026 Catalogue





“ The Covetrus Branded Product portfolio continues to grow. From high-performing, in-practice supplies to general health and wellness products, we’ve got you covered. ”



Covetrus Branded Products

Covetrus is committed to offering a reliable range of products, enabling professionals to dedicate their attention to patient care, medical practice, and caregiving. Covetrus Branded Products are designed to deliver quality and value at competitive prices, supporting both small and large animal practices in enhancing their operational efficiency, profitability and clinical care.

Our comprehensive portfolio features Covetrus Essentials and our retail line, Pro Pet Solutions. Covetrus Essentials provides clinics and practices with high-quality, everyday essentials at exceptional value, while Covetrus Pro Pet Solutions offers pet owners veterinarian-recommended, affordable products for optimal pet care.

We invite you to explore the breadth and potential of our product assortment, as we continue to expand our offerings across all categories to meet the evolving needs of veterinary professionals and pet owners.

The products presented in this catalogue are intended for veterinary use and in animals under veterinary guidance. For use in animals only. They are not intended or suitable for human use.

Covetrus Essentials	2
1 Injection	7
Syringes	7
Reinforced Needles	8
Syringe and Needle Combination	8
Hypodermic Needles	9
Supplementary	10
2 Infusion & Transfusion	11
IV Catheters	11
Obturator	11
Winged Needles	12
IV Infusion Sets	12
Infusion Set Accessories	12
Miscellaneous	12
3 Surgery & Anaesthesia	13
Two-layer Non-woven Surgical Drapes	13
One-layer PE Surgical Drapes	14
Textile Surgical Drapes	14
Vet Surgical Drape Sets	14
Accessories	15
Autoclave	15
Sterilisation	15
Anaesthesia	15
4 Cleaning & Disinfection	17
Hand & Skin	17
Surface Disinfection	18
Instrument Disinfection	20
Instrument Cleaning & Maintenance	21
Complementary Products	21
Hygiene Paper	22
5 Personal Protective Equipment (PPE)	23
Clothes	23
Shoe Covering	23
X-ray Protective Clothing	24
Latex Examination Gloves	24
Nitrile Examination Gloves	25
Vinyl Gloves	25
Veterinary Gloves	26
Surgical Gloves	26
Face Masks	26
6 Sutures	27
Suture Summary	27
Spool Sutures	32
Suture Needles	35
Surgical Skin Adhesive	38
Skin Staplers	38
7 Wound Care, Skin Care & Bandages	39
Cohesive Bandage	39
Orthopaedic Casting	40
Bandage Covering	40
Tape	41
Padding	42
Other Bandages	43
Non-woven	44
Hydrophilic	44

Gauze & Cotton Balls	45
Wound Care	45
Skin Care	46
Miscellaneous	46
Hot/Cold	46
8 Pharma & Chemicals	47
Prescription Vials	47
Pharmaceuticals	47
Chemicals - Water	50
Chemicals - Alcohol	51
9 Urology	51
Urinary Catheters	51
Urine Collection	52
Urine Container	52
10 Diagnostics	53
Urine Test Strips	53
Faeces Analysis	53
Dermatophyte Test	54
Rapid Tests	55
Physical Examination	56
Thermometers	57
Percussion Hammer	57
Penlights	57
Biopsy Punch	58
Dermal Curette	58
Miscellaneous	58
11 Laboratory	59
Microscope Slides	59
Pipette	59
Applicator	59
Microscope Cleaning Set	59
12 Patient Care & Animal Handling	60
Identification	60
Microchip Reader	60
Protective Gear	61
Muzzles	61
Ectoparasites	62
Nail Clippers	63
Miscellaneous	63
Underpads	63
Cadaver Bags	63
Cat Litter	64
13 Large Animals	65
Cow Ruminal Magnets	65
Marking	65
Veterinary Gels	66
Teat Instruments	66
14 Equipment	67
Patient Monitoring/Surgery	67
Handheld Veterinary Monitor	68
Multi-Parameter Patient Monitor	69
Miscellaneous	71
Veterinary Glucose Meter	71
Handheld Refractometer	71
Centrifuges	72
Digital Infusion Devices	73
Vet Piezo Scaler	74

	Handpieces	75
	Straight Handpieces 1:1	75
	Water Attachments	76
	LED Curing Light	76
	Thermosealer	76
	Ultrasonic Cleaner	76
	Microscopes	77
	Scales	78
	Clippers and Trimmers	78
15	Instruments	79
	Needle Holders	79
	Premium Range with TC Insert	80
	Scalpel Handles & Blades	81
	Scalpel Blades	82
	Scissors	83
	Forceps	87
	Towel Clamps	95
	Wound Retractors	96
	Bone Instruments	96
	Other Products	97
	Instruments for Large Animals	98
	Instrument Boxes & Accessories	98
	Maintenance	100
16	Dentistry	101
	Hygiene & Maintenance	101
	Dental Care	101
	Dental Instruments	102
	Parodontology	107
	Application & Modelling Tools	108

	Covetrus Pro Pet Solutions	109
1	Nutraceuticals	110
	Joint support	110
	Urinary support	111
	Skin & Coat support	111
	Immune support	112
	Liver support	113
	Stress support	113
	Puppy support	114
	Hairball Remedy	114
	Kidney support	114
	Intestinal support	115
	Rehydration support	115
	Gastro support	116
	Metabolism support	116
	Dental support	116
	Other dental chew	116
2	Topicals and Shampoo bars	117
	Topicals	117
	Shampoo bars	118
3	Patient Care & Animal Handling	119
	Protective boots	119
	Cat Litter	119
	Muzzles	120
	Ectoparasites	120

VISIT
products.covetrus.com
TO LEARN MORE



Covetrus Essentials Best Sellers

MOST WANTED | ESSENTIALS

Personal Protective Equipment (PPE)

Nitrile Powder-free Gloves

- Powder-free ensuring an optimal barrier protection
- Finger textured surface
- Different colours & types available
- Ambidextrous
- Beaded cuff



PAGE 25



Injection

Disposable Syringes

- Made of the highest quality plastic and have excellent tolerances for the perfect extrusion of liquids.

PAGE 7



Equipment

Pet Hair Clipper

- Compact design
- Low noise
- Rounded Blade Teeth
- Detachable ceramic Blades
- Strong cutting power
- Waterproof & USB Rechargeable

PAGE 78





Covetrus Brand Most Wanted Products

A selection of our best products across a variety of in-clinic supply and veterinary categories.

MOST WANTED | PRO PET SOLUTION

Hairball Remedy

CatLax Hairball Remover

- CatLax contains ingredients that gently lubricate the digestive tract to prevent hairball build-up
- The paste acts as a lubricant that coats the swallowed hair and allows the hairball to travel down the digestive tract
- It has a mild laxative effect that penetrates hardened stool so it can be passed



PAGE 114

Topicals

Eye Drops

- For dry eyes, redness, and irritation
- Preservative-free Eye Drops
- Soothes and refreshes
- 100% natural: with chamomile, euphrasia and hamamelis
- Content: 10ml



PAGE 117

Shampoo bars

Pet Shampoo Bars

- Made from natural and naturally-derived ingredients
- pH-balanced for a clean and gentle wash
- No plastic packaging
- 100% vegan
- Eco friendly
- Degradable without residue

PAGE 118




Injection

The Covetrus Syringe - Portfolio of Products

Syringes

Covetrus Syringes are made of the highest quality plastic and have excellent tolerances for the optimal extrusion of liquids. The scale on the syringe makes an optimal dosage possible.

2-part Syringes
The 2 ml has a scale up to 3 ml, the 5 ml up to 6 ml, the 10 ml up to 12 ml, the 20 ml up to 24 ml.



- Luer connection
- Disposable, sterile

Item	Description	Quantity
1	2802508 Syringe Disposable 2-part 1 ml, Luer Centric	100 pcs
2	2802509 Syringe Disposable 2-part 2 ml, Luer Centric	100 pcs
3	2802510 Syringe Disposable 2-part 5 ml, Luer Eccentric	100 pcs
4	2802511 Syringe Disposable 2-part 10 ml, Luer Eccentric	100 pcs
5	2802512 Syringe Disposable 2-part 20 ml, Luer Eccentric	50 pcs


3-part Syringes
• Luer connection
• Available in Centre or Side (eccentric) position of the conus



Item	Description	Quantity
1	2802548 Syringe Disposable 3-part 1 ml, Luer Centric with Low Dead Space Plunger	100 pcs
2	2802518 Syringe Disposable 3-part 1 ml, Luer Centric	100 pcs
3	2802513 Syringe Disposable 3-part 3 ml, Luer Centric	100 pcs
4	2802514 Syringe Disposable 3-part 5 ml, Luer Centric	100 pcs

Item	Description	Quantity
1	2802519 Syringe Disposable 3-part 5 ml, Luer Eccentric	100 pcs
2	2802515 Syringe Disposable 3-part 10 ml, Luer Eccentric	100 pcs
3	2802516 Syringe Disposable 3-part 20 ml, Luer Eccentric	50 pcs
4	2802521 Syringe Disposable 3-part 30 ml, Luer Eccentric	50 pcs
5	2802517 Syringe Disposable 3-part 50 ml, Luer Eccentric	50 pcs

Injection



3-part Syringes

- Luer Lock connection
- Centre position of the conus

Item	Description	Quantity
1	2802522 Syringe Disposable 3-part 3 ml, Luer Lock Centric	100 pcs
2	2802523 Syringe Disposable 3-part 5 ml, Luer Lock Centric	100 pcs
3	2802524 Syringe Disposable 3-part 10 ml, Luer Lock Centric	100 pcs



Insulin Syringes

3-part Insulin Syringes with integrated needle. Made of the highest quality plastic and have excellent tolerances for the optimal extrusion of liquids. Sterile, disposable.

Item	Description	Quantity
1	2802549 Insulin Syringe Disposable U100 30G x 1/2", (0.3 x 12 mm) 1ml	100 pcs
2	2802550 Insulin Syringe Disposable U40 30G x 1/2", (0.3 x 12 mm) 1ml	100 pcs

Reinforced Needles



High quality reinforced disposable injection needle. Specially designed for the veterinary profession. Thickened needle hub for extra strength and to avoid break of needle hub. Made of high-quality stainless steel and polypropylene. Sterilised by ethylene oxide.


- Luer lock
- Non-pyrogenic
- The needle hub strength has been tested at 100-110N
- Made of high-quality stainless steel and polypropylene

Item	Description	Dimensions (G x inch)	Dimensions (mm)	Colour	Quantity
2801486	Reinforced Disposable Needles	17G x 2"	1.5 x 50	Pink	100 pcs
2801489	Reinforced Disposable Needles	15G x 1 1/2"	1.8 x 38	Dark grey	100 pcs
2801487	Reinforced Disposable Needles	14G x 1 1/2"	2.1 x 38	Teal	100 pcs
2801488	Reinforced Disposable Needles	14G x 2 1/2"	2.1 x 60	Teal	100 pcs
2801947	Reinforced Disposable Needles	14G x 3 1/4"	2.1 x 80	Teal	100 pcs

Needle and Syringe Combination

Disposable syringe with fixed needle. 3-component syringe with centric luer adaptor. Delicately designed, facet shaped needle point with features of sharpness, smooth puncture, less damage to the tissue, less painful feeling to the patient.

- Semi-transparent plastic needle hub for the convenient
- Observation of blood reflux



Extra sharp facet shaped needles



Item	Description	Syringe	Dimensions (G x inch)	Dimensions (mm)	Colour	Quantity
2801951	Disposable Syringe/Needle	2-3 ml	23G x 1"	0.6 x 25	Blue	100 pcs
2801950	Disposable Syringe/Needle	2-3 ml	21G x 5/8"	0.8 x 16	Green	100 pcs
2801952	Disposable Syringe/Needle	5-6 ml	21G x 1 1/2"	0.8 x 40	Green	100 pcs
2801949	Disposable Syringe/Needle	5-6 ml	18G x 1 1/2"	1.25 x 40	Pink	100 pcs
2801948	Disposable Syringe/Needle	20-25 ml	18G x 1 1/2"	1.25 x 40	Pink	50 pcs


Hypodermic Needles

Hypodermic needles for single use. The hypodermic needles stand for optimised siliconisation of needle tube, high pull-out forces of cannula-hub-connection, state of the art anti-coring treatment of needle tube, precision lancet point according to ISO 9626 with exact smooth bevel cut and consistent low penetration load of needle and skin.

- A range of sizes specially designed for the veterinary profession
- Siliconised needle tubes for improved slippage
- Made of high quality stainless steel
- Delicately designed, facet shaped needle point with features of sharpness, smooth puncture, less damage to the tissue, less painful feeling to the patient
- Semi-transparent plastic needle hub for the convenient observation of blood reflux with ISO colour-code for quick identification of needle diameter
- Sterilised by ethylene oxide



A short bevel type has a greater angle compared to a regular bevel and because of this, it may facilitate the injection of fluids with high-viscosity.



Extra sharp facet shaped needles

Item	Description	Dimensions (G x inch)	Dimensions (mm)	Colour	Quantity
2802547	Hypodermic Needle Long Bevel Luer Lock	30G x 1/2"	0.30 x 12	Yellow	100 pcs
2802546	Hypodermic Needle Long Bevel Luer Lock	27G x 3/4"	0.40 x 20	Grey	100 pcs
2802545	Hypodermic Needle Long Bevel Luer Lock	26G x 1/2"	0.45 x 12	Brown	100 pcs
2802544	Hypodermic Needle Long Bevel Luer Lock	26G x 1"	0.45 x 25	Brown	100 pcs
2802543	Hypodermic Needle Long Bevel Luer Lock	25G x 5/8"	0.50 x 16	Orange	100 pcs
2802542	Hypodermic Needle Long Bevel Luer Lock	24G x 1"	0.55 x 25	Violet	100 pcs
2802541	Hypodermic Needle Long Bevel Luer Lock	23G x 1"	0.60 x 25	Blue	100 pcs
2802540	Hypodermic Needle Long Bevel Luer Lock	23G x 1 1/4"	0.60 x 30	Blue	100 pcs
2802539	Hypodermic Needle Long Bevel Luer Lock	23G x 5/8"	0.60 x 16	Blue	100 pcs
2802538	Hypodermic Needle Long Bevel Luer Lock	23G x 3/4"	0.60 x 20	Blue	100 pcs
2802536	Hypodermic Needle Long Bevel Luer Lock	22G x 3/4"	0.70 x 20	Black	100 pcs
2802537	Hypodermic Needle Long Bevel Luer Lock	22G x 1 1/4"	0.70 x 30	Black	100 pcs
2802533	Hypodermic Needle Long Bevel Luer Lock	21G x 5/8"	0.80 x 16	Green	100 pcs
2802535	Hypodermic Needle Long Bevel Luer Lock	21G x 1"	0.80 x 25	Green	100 pcs
2802534	Hypodermic Needle Long Bevel Luer Lock	21G x 1 1/2"	0.80 x 40	Green	100 pcs
2802532	Hypodermic Needle Long Bevel Luer Lock	20G x 1"	0.90 x 25	Yellow	100 pcs
2802531	Hypodermic Needle Long Bevel Luer Lock	20G x 1 1/2"	0.90 x 40	Yellow	100 pcs
2802530	Hypodermic Needle Long Bevel Luer Lock	19G x 1"	1.10 x 25	Beige	100 pcs
2802528	Hypodermic Needle Long Bevel Luer Lock	18G x 1"	1.20 x 25	Pink	100 pcs
2802529	Hypodermic Needle Long Bevel Luer Lock	18G x 1 1/2"	1.20 x 40	Pink	100 pcs
2802527	Hypodermic Needle Short Bevel Luer Lock	18G x 1 1/2"	1.20 x 40	Pink	100 pcs
2802525	Hypodermic Needle Short Bevel Luer Lock	16G x 1"	1.60 x 25	White	100 pcs
2802526	Hypodermic Needle Short Bevel Luer Lock	16G x 1 1/2"	1.60 x 40	White	100 pcs

Supplementary

Withdrawal Cannula Spike Plus

Withdrawal cannula with snap fit cap for withdrawing fluids from vials or semi-rigid containers. Slim, short spike, for effortless piercing of even thicker stoppers.

- Reduces contamination
- Available with green colour caps. PVC-free, DEHP-free
- Luer Lock



Item	Description	Quantity
2800933	Withdrawal Cannula Spike Plus Green	100 pcs

Syringe Stopper

Syringe stopper from medical silicone.



Item	Description	Quantity
2802562	Syringe Stopper, Non-sterile	100 pcs

Sharps Container

Single use container for safe disposal of clinical waste. Ability to permanently lock for final disposal.

- Marking for maximum filling level
- Autoclavable
- The lid can be sealed either temporarily or permanently
- Manufactured from environmentally friendly PP material
- Hypodermic needle remover integrated in the lid

Item	Description	Size
1	2800203	Sharps Container 0.25 L
2	2800204	Sharps Container 0.5 L
3	2800205	Sharps Container 1L
4	2800206	Sharps Container 1.5 L
5	2800207	Sharps Container 2L
6	2800208	Sharps Container 2.2 L
7	2800209	Sharps Container 3L
8	2800210	Sharps Container 3.5 L
9	2800211	Sharps Container 4 L
10	2800212	Sharps Container 6 L
11	2800213	Sharps Container 10 L



Infusion & Transfusion

IV Catheters

IV Catheter with Port & Wings



Item	Colour	Dimensions (G x inch)	Dimensions (mm)	Quantity
2800707	Purple	26G x 3/4"	0.64 x 19 mm	50 pcs
2800706	Yellow	24G x 3/4"	0.7 x 19 mm	50 pcs
2800705	Blue	22G x 1"	0.9 x 25 mm	50 pcs
2800704	Pink	20G x 1 1/4"	1.0 x 32 mm	50 pcs

IV Catheter without Port & with Wings



Item	Colour	Dimensions (G x inch)	Dimensions (mm)	Quantity
2800721	Yellow	24G x 3/4"	0.75 x 19 mm	50 pcs
2800720	Blue	22G x 1"	0.9 x 25 mm	50 pcs
2800719	Pink	20G x 1 1/4"	1.0 x 32 mm	50 pcs
2800718	Green	18G x 1 3/4"	1.3 x 45 mm	50 pcs
2800717	Dark Green	18G 1 1/4"	1.3 x 32 mm	50 pcs
2800716	Grey	16G x 1 3/4"	1.75 x 45 mm	50 pcs
2800715	Orange	14G x 1 3/4"	2.0 x 45 mm	50 pcs

IV Catheter for Large Animals without Port and without Wings



Item	Dimensions (G x inch)	Dimensions (mm)	Quantity
2800722	14G x 2 3/4"	2.00 x 70 mm	30 pcs

IV Catheter without Port and without Wings



Item	Colour	Dimensions (G x inch)	Dimensions (mm)	Quantity
2800714	Yellow	24G x 3/4"	0.75 x 19 mm	50 pcs
2800713	Blue	22G x 1"	0.90 x 25 mm	50 pcs
2800712	Pink	20G x 1 1/4"	1.00 x 32 mm	50 pcs

Obturators

Obturators are recommended to avoid the formation of clots in the cannula lumen after each suspension of the IV therapy.

Obturator for Dual Entry IV Catheter

Obturators are recommended to avoid the formation of clots in the cannula lumen after each suspension of the IV therapy.



Item	Colour	Dimensions	Quantity
2800723	Yellow	24G	50 pcs
2800724	Blue	22G	50 pcs
2800725	Pink	20G	50 pcs

Obturator for Single Entry IV Catheter



Item	Colour	Dimensions	Quantity
2800726	Yellow	24G	50 pcs
2800727	Blue	22G	50 pcs
2800728	Pink	20G	50 pcs

Winged Needles

Winged Needles Infusion Set

Winged infusion set also known as “butterfly” or “scalp vein set”.



Item	Colour	Dimensions (G x inch)	Dimensions (mm)
2800693	●	27G x 3/4"	0.4 x 20 mm
2800694	●	25G x 3/4"	0.5 x 20 mm
2800695	●	23G x 3/4"	0.6 x 20 mm
2800696	●	22G x 3/4"	0.7 x 20 mm
2800697	●	21G x 3/4"	0.8 x 20 mm
2800698	●	19G x 3/4"	1.1 x 20 mm
2800699	●	18G x 3/4"	1.2 x 20 mm

IV Infusion Sets

IV Infusion Set with Rotating Luer Lock Connection

IV Infusion set with airvent and Macrodrop 20 drops/ml. Rotating luer lock or luer slip connection. Including an Y injection site or latex injection bulb for injecting drugs in the IV system.



Item	Length	Injection site	Filter	Drops	Quantity
2801866	180 cm	Injection Bulb	15 µm	20 ml	30 pcs
2801867	180 cm	No Injection Bulb	15 µm	20 ml	30 pcs
2801868	150 cm	Injection Bulb	No Filter	20 ml	30 pcs
2801869	150 cm	Injection Bulb	No Filter	20 ml	30 pcs
2801870	200 cm	Y-injection Site	15 µm	20 ml	30 pcs
2801871	200 cm	Y-injection Site	No Filter	20 ml	30 pcs

Infusion Set Accessories

Heidelberg extension line

Disposable extension line, low pressure, non-toxic, non-pyrogenic, single use.



Item	Dimensions	Quantity
2801872	30 cm	100 pcs
2801873	75 cm	100 pcs
2801874	100 cm	50 pcs
2801875	140 cm	50 pcs

Combi Stopper

Universal cap for closing luer and luer lock connections. Suitable for closure of IV cannulas, IV sets, and central venous catheters.



Item	Colour	Quantity
2800937	●	300 pcs

3-Way Stop Cock

Usable with multiple infusions, featuring a 360° rotatable male luer and two female luer connectors.



Item	Colour	Quantity
2801829	●	50 pcs
2801830	●	50 pcs

Miscellaneous

Tourniquet Small Animals

Extra slim tourniquet, ideal for cats and small dogs, latex-free.



Item	Colour	Quantity
2801885	●	1 pc

Tourniquet Regular

Tourniquet with patented closing mechanism that can be operated with one hand, latex-free.



Item	Colour	Quantity
2801011	●	1 pc



Surgery & Anaesthesia

Two-layer Non-woven Surgical Drapes

Adhesive Fenestrated Drape

- Sterile, individually wrapped
- Adhesive oval fenestration
- Two layer material: polyethylene film with absorbent non-woven top layer
- Impermeable to moisture and micro-organisms



Item	Dimensions	Quantity
2800085	30 x 45 cm, Ø 2 x 6.5 cm	1 pc
2800086	30 x 45 cm, Ø 5 x 10 cm	1 pc
2800087	30 x 45 cm, Ø 7.5 x 15 cm	1 pc
2800091	45 x 75 cm, Ø 2 x 6.5 cm	1 pc
2800088	45 x 75 cm, Ø 9 x 20 cm	1 pc
2800089	45 x 75 cm, Ø 8 x 25 cm	1 pc
2800090	88 x 75 cm, Ø 27 x 13 cm	1 pc

Adhesive Fenestrated Drape 2-part

- Sterile, individually wrapped
- Variable adhesive fenestration, two parts can be put on top of each other to create a smaller fenestration
- Two layer material: polyethylene film with absorbent non-woven top layer
- Impermeable to moisture and micro-organisms



Item	Dimensions	Quantity
2800101	45 x 75 cm, Ø 7 cm	1 pc

Non-adhesive Fenestrated Drape

- Sterile, individually wrapped
- Two layer material: polyethylene film with absorbent non-woven top layer
- Impermeable to moisture and micro-organisms

Item	Dimensions	Quantity
2800099	50 x 70 cm, Ø 7 cm	1 pc
2800100	75 x 90 cm, Ø 10 cm	1 pc

Adhesive Drape Non-fenestrated

- Sterile, individually wrapped
- With 5 cm adhesive strip
- Two layer material: polyethylene film with absorbent non-woven top layer
- Impermeable to moisture and micro-organisms



Item	Dimensions	Quantity
2800094	50 x 50 cm + 5 cm adhesive strip	1 pc
2800092	50 x 70 cm + 5 cm adhesive strip	1 pc
2800093	75 x 90 cm + 5 cm adhesive strip	1 pc

Non-adhesive Non-fenestrated Drape

- Sterile, individually wrapped
- Two layer material: polyethylene film with absorbent non-woven top layer
- Impermeable to moisture and micro-organisms




Item	Dimensions	Quantity
2800095	40 x 50 cm	1 pc
2800096	50 x 50 cm	1 pc
2800097	50 x 70 cm	1 pc
2800098	75 x 90 cm	1 pc



Textile Surgical Drapes

Textile Surgical Drapes

- 100% cotton, non sterile



Item	Dimensions	Quantity
2800675	36 x 60 cm	1 pc
2800676	50 x 50 cm	1 pc
2800677	70 x 70 cm	1 pc
2800678	90 x 90 cm	1 pc
2800679	100 x 100 cm	1 pc
2800680	100 x 150 cm	1 pc
2800681	200 x 250 cm	1 pc

Textile Drapes with Slit

- 100% cotton, non sterile



Item	Dimensions	Quantity
2800684	40 x 60 cm, slit (15 cm)	1 pc
2800685	50 x 50 cm, slit (10.5 cm)	1 pc
2800686	70 x 70 cm, slit (20 cm)	1 pc
2800687	90 x 90 cm, slit (35 cm)	1 pc
2800688	100 x 100 cm, slit (20 cm)	1 pc
2800689	100 x 150 cm, slit (22 cm)	1 pc
2800690	200 x 250 cm, slit (25 cm)	1 pc

Textile Fenestrated Drapes

- 100% cotton, non sterile



Item	Dimensions	Quantity
2800682	∅ 7 cm, 50 x 50 cm	1 pc
2800683	∅ 7.5 cm, 100 x 100 cm	1 pc

One-layer Non-woven Surgical Drapes

One-layer PE Surgical Drape

Made of soft one-layer polyethylene.

- Tear resistant
- Impermeable to moisture and micro-organisms
- Sterile



Item	Dimensions	Quantity
2800794	30 x 45 cm	1 pc
2800795	45 x 60 cm	1 pc
2800796	60 x 90 cm	1 pc
2800797	90 x 120 cm	1 pc
2800798	120 x 120 cm	1 pc
2800799	120 x 180 cm	1 pc
2800800	120 x 250 cm	1 pc

Vet Surgical Drape Sets

Set Small Animals

The sterile vet surgery set can be used for procedures on small animals. The drapes are made of soft two-layer material: polyethylene film with absorbent non woven top layer.

- The set contains:
- 1 x Mayo cover 80 x 145 cm
 - 4 x Surgical drape 75 x 75 cm with adhesive strip
 - 20 x Compress
 - 1 x Surgical drape 200 x 240 cm
 - 1 x Table cover 150 x 190 cm
 - 1 x Wrapping sheet 80 x 80 cm



Item	Quantity
2800103	1 pack

Set Small/Dental Surgery

The sterile small surgery set contains drapes for oral surgery and implantology procedures.

- The set contains:
- 2 x Surgical drape 45 x 75 cm
 - 1 x Table cover 75 x 90 cm
 - 1 x Dental drape 160 x 150 cm
 - 2 x Cord drape for tubing 8 x 120 cm
 - 2 x OP adhesive tape 10 x 50 cm
 - 1 x Wrapping sheet 75 x 80 cm



Item	Quantity
2800079	1 pack

Colic Surgery Set

The sterile horse colic surgery set contains drapes to cover the horses body and the hind legs with 2 bags attached to hold the intestines. The sterile drapes are made of soft two-layer material: polyethylene film with absorbent non-woven (spunbond) top layer.

- The set contains:
- 2 x Stockinette 75 x 120 cm
 - 1 x Perforated drape 260 x 370 cm
 - 1 x Table cover 150 x 200 cm



Item	Quantity
2800080	1 pack

Accessories

Adhesive Operation Tape

Sterile blue tape with adhesive to keep the drape in place.



Item	Dimensions	Quantity
2800102	10 x 50 cm	1 pc

Autoclave

Autoclave Tape

Self-adhesive autoclave tape for steam and chemical vapour sterilisers. Changes colour to show sterilisation has taken place.



Item	Dimensions	Quantity
2800887	12 mm x 50 m	1 pc
2800888	19 mm x 50 m	1 pc
2800889	25 mm x 50 m	1 pc

Autoclave Indicator Marker

Marker with chemical indicator ink for steam sterilization. Marks on labels or pouch bags change color from wine-red to green after sterilization.



Item	Quantity
2800860	1 pc

Sterilisation

Sterilisation Roll

Sterilisation roll for making sterilisation pouches using a sealing device. With triple sealing seam and indicator field. Sterilisation rolls maintain sterility of device for up to 180 days after sterilisation, if stored correctly.



Item	Dimensions	Quantity
2802345	50 mm x 200 m	1 roll
2802346	75 mm x 200 m	1 roll
2802347	100 mm x 200 m	1 roll
2802348	150 mm x 200 m	1 roll
2802349	200 mm x 200 m	1 roll
2802350	250 mm x 200 m	1 roll
2802351	300 mm x 200 m	1 roll

Self-sealing sterilization pouch

Self-seal sterilisation pouches made from polyethylene and paper. With triple seal seam, wide adhesive strip and indicator field. Suitable for autoclaving or chemical sterilisation. Can be heated to 135 degrees Celsius. Colour indicator changes blue to brown when sterilisation has taken place.



Item	Dimensions	Quantity
2802340	83 x 159 mm	200 pcs
2802341	89 x 254 mm	200 pcs
2802342	134 x 280 mm	200 pcs
2802343	190 x 356 mm	200 pcs
2802344	254 x 381 mm	200 pcs

Anaesthesia

Anaesthesia Gas Filter Canister

Adsorbs up to 300g of anaesthetic agent.



Item	Quantity
2802505	1200 g

Neoprene Re-Breathing Bag

Neoprene re-breathing bags are high-quality, latex-free reservoirs designed for veterinary anaesthesia systems. They ensure reliable gas delivery and allow clear visualization of respiratory movements.



Item	Quantity
2802626	0,5 L
2802627	1 L
2802628	2 L
2802629	3 L

Latex Re-Breathing Bag

The latex re-breathing bags are flexible, durable reservoirs. They provide controlled gas flow and enable easy monitoring of respiratory effort. Additionally, they come with a hanging loop.



Item	Quantity
2802630	0,5 L
2802631	1 L
2802632	2 L
2802633	3 L

Veterinary Anaesthesia Machine CA23

The universal veterinary anaesthesia machine is applicable to clinical anaesthesia of dogs, cats, pigs, monkeys, rodent, reptiles, birds etc. With minimal operation design, vets shall feel easy, safe and reliable to perform surgery on pets.



Item	Description
2802314	Veterinary Anaesthesia Machine

Standard Configuration (Scope of delivery)

Item	Description
2802499	Vaporizer for Isoflurane - Key Fill - Cagemount
2802500	Vaporizer Endcaps (1 pair)
2802494	Breathing Circuit assembly with APL valve + pressure gauge
2802491	Breathing Circuit Duo-limb, 22mm(M)/15mm(F)
2802492	Mapleson F Non-rebreathing Circuit (Jackson-Rees)
2802493	Corrugated Tube - ID 22.0mm
2802496	Air Hose(3-layer Compounded Rubber) - ID8.0-9/16(SGL)
2802497	Custom-made-Air Hose (6 mm) 2,0 m
2802498	Custom-made-tube holder
2802495	CO2 absorber, (size 153*99*148mm) 1.100ml

Optional Configuration

Item	Description
2802488	Wall mounting bracket
2802501	Wall Hanging Frame
2802489	R-AB H(H-base)
2802490	R-HS(Oxygen Concentrator Shelf for H-base)

Endotracheal Tube

Disposable endotracheal tube with low pressure cuff and standard connector. Sterile, Murphy type.

Item	I.D.	CH	Length	Quantity
2800752	2 mm	10	14.5 cm	10 pcs
2800753	2.5 mm	12	15 cm	10 pcs
2800754	3 mm	14	16.5 cm	10 pcs
2800755	3.5 mm	16	18.5 cm	10 pcs
2800756	4 mm	18	21 cm	10 pcs
2800757	4.5 mm	18	22.5 cm	10 pcs
2800758	5 mm	20	24.5 cm	10 pcs
2800759	5.5 mm	22	27.5 cm	10 pcs
2800760	6 mm	24	28.5 cm	10 pcs
2800761	6.5 mm	26	29.5 cm	10 pcs
2800762	7 mm	28	31 cm	10 pcs
2800763	7.5 mm	30	31.5 cm	10 pcs
2800764	8 mm	32	33 cm	10 pcs
2800765	8.5 mm	34	33 cm	10 pcs
2800766	9 mm	36	33 cm	10 pcs
2800767	9.5 mm	38	33 cm	10 pcs
2800768	10 mm	40	33 cm	10 pcs
2800769	10.5 mm	42	33 cm	10 pcs
2800770	11 mm	44	33 cm	10 pcs



Cleaning & Disinfection

Hand and Skin

Hand Lotion

- Delivers long-lasting, noticeable effect even with a small amount
- Provides replenishing and moisturizing benefits

pH value 6.5



Item	Packaging type	Quantity
2802352	Tube	100 ml

Hand Wash Lotion

- Prevents hands from drying out, even with frequent washing
- Cleanses and nourishes

pH value 6.5



Item	Packaging type	Quantity
2802353	Bottle	1L
2802354	Canister	5L

Skin and Hand Disinfection

- Excellent skin compatibility, dermatologically tested as "very good"
- Short exposure time with broad-spectrum efficacy on polyalcohol basis
- Broad-spectrum efficacy: bactericidal, levurocidal, tuberculocidal, and limited virucidal

pH value 5.5



Item	Packaging type	Quantity
2802355	Bottle	150 ml
2802356	Bottle	1L

Hand Disinfection Gel

- "3-in-1" formulation: non-drip gel, liquefies on contact for easy spread, and leaves lasting skin care
- Dermatologically tested with "very good" rating
- Extremely short exposure time with broad-spectrum efficacy (bactericidal, levurocidal, tuberculocidal, limited virucidal)

pH value 5.5



Item	Packaging type	Quantity
2802357	Dispenser	250 ml
2802358	Bottle	1L

Hand Wash Lotion

- Gentle skin wash with balanced composition and pH identical to skin, suitable for all skin types in hospital settings
- Free from aldehydes and phenols, ensuring safe and non-irritating use



Item	Packaging type	Quantity
2802461	Bottle	500 ml

Alco Aloe Hydroalcoholic Gel

- Ready-to-use antiseptic gel for hand disinfection
- Broad-spectrum efficacy: Bactericidal and fungicidal action for reliable protection



Item	Packaging type	Quantity
2802465	Bottle	500 ml



Some products from this chapter may be available only in selected countries due to the local registration/regulatory requirements. Please check availability with your Covetrus representative.

Surface Disinfection

Surface Disinfection Wipes

- Excellent material compatibility
- Highly effective rapid disinfection with significantly reduced alcohol content
- 13 x 20 cm PP/VIS mixed fabric wipes (120 pieces), ideal for small and hard-to-reach surfaces
- Broad-spectrum efficacy: bactericidal, levurocidal, tuberculocidal, and limited virucidal PLUS



Item	Packaging type	Quantity
2802359	Refill pack	120 pcs
2802360	Container	120 pcs

Surface Disinfection Dry Wipes

- High-quality, thick, lint-free fleece material with excellent liquid absorption and release
- Economical in use and pleasantly soft to the touch
- Size: 28 x 30 cm, 90 wipes, made of 100% PET



Item	Packaging type	Quantity
2802361	Roll	90 pcs
2802478	Dispenser	1

Surface Disinfection

- Ready-to-use rapid disinfectant for small surfaces with excellent material compatibility (e.g., tablets, acrylic glass, delicate synthetic leather)
- Strongly reduced alcohol content while maintaining top-level efficacy
- Broad-spectrum efficacy: bactericidal, levurocidal, tuberculocidal, and limited virucidal PLUS



Item	Packaging type	Quantity
2802362	Bottle	1L
2802363	Canister	5L

Surface Disinfection FV

- Ready-to-use, alcohol-free solution with innovative formulation
- Fast-acting and virucidal
- Also suitable as soaking liquid for wet wipes
- Ready-to-use, alcohol-free solution with innovative formulation
- Broad-spectrum efficacy: bactericidal, levurocidal, tuberculocidal, and virucidal



Item	Packaging type	Quantity
2802364	Bottle	1L
2802365	Canister	5L

Surface Disinfection Spray, non-alcoholic

- Ready-to-use aqueous solution for disinfecting
- Excellent material compatibility – even for alcohol-sensitive surfaces and materials
- Broad-spectrum efficacy: bactericidal, levurocidal, tuberculocidal, and limited virucidal PLUS



Item	Packaging type	Quantity
2802366	Bottle	1L
2802367	Canister	5L

Surface Disinfection Concentrate, non-alcoholic

- Aldehyde-free concentrate for surface disinfection and cleaning
- Excellent cleaning performance without additional additives
- Leaves no sticky residues
- Broad-spectrum efficacy: bactericidal, levurocidal, tuberculocidal, and limited virucidal PLUS



Item	Packaging type	Quantity
2802368	Canister	5L

Surface Disinfection Wipes

- Excellent material compatibility, ideal for sensitive environments
- Aldehyde- and phenol-free formulation
- Broad-spectrum efficacy: Bactericidal, yeasticidal, and general virucidal action for comprehensive hygiene



Item	Packaging type	Quantity
2802435	Container	100 pcs

Surface Disinfection Wipes, non-alcoholic

- Alcohol-free, ready-to-use wipes for rapid cleaning and disinfection of medical surfaces and devices
- Formulated with hydrogen peroxide and quaternary ammonium, ensuring material safety without fading or staining
- Broad-spectrum efficacy: bactericidal, fungicidal, virucidal, and sporicidal



Item	Packaging type	Quantity
2802467	Container	100 pcs

Surface Disinfection

Alcohol-based, ready-to-use disinfectant for rapid cleaning of surfaces

- Free from aldehydes and phenols, safe for materials without fading or staining
- Broad-spectrum efficacy: bactericidal, fungicidal, and virucidal



Item	Packaging type	Quantity
2802469	Bottle	1L
2802470	Canister	5L

Hypochlorous Surface Disinfection

- Alcohol-free surface disinfectant with excellent material compatibility
- Broad-spectrum efficacy: bactericidal, yeasticidal, and general virucidal



Item	Packaging type	Quantity
2802471	Bottle	1L
2802472	Canister	5L

Surface Disinfection, non-alcoholic

- Alcohol-free, ready-to-use disinfectant for rapid cleaning and disinfection of surfaces
- Formulated with hydrogen peroxide and quaternary ammonium, safe for materials without discoloration or staining
- Broad-spectrum efficacy: bactericidal, fungicidal, virucidal, and sporicidal



Item	Packaging type	Quantity
2802436	Bottle	1L
2802468	Canister	5L



Some products from this chapter may be available only in selected countries due to the local registration/regulatory requirements. Please check availability with your Covetrus representative.

Instrument Disinfection

Enzymatic Instrument Cleaner, Concentrate

- Low-foam enzymatic cleaner for medical instruments, endoscopes, and lab equipment
- 4-enzyme formula: Protease, Amylase, Lipase, Cellulase
- trasonic bath compatible



Item	Packaging type	Quantity
2802369	Bottle	1L

Instrument Disinfection, Concentrate FV

- Liquid concentrate for manual disinfection of medical and surgical instruments using the immersion bath method
- Free from peracids, non-corrosive yet highly effective, ensuring excellent material compatibility
- Broad-spectrum action: bactericidal, tuberculocidal, fungicidal, virucidal



Item	Packaging type	Quantity
2802370	Bottle	2L

Enzymatic Instrument Cleaner, Concentrate

- Ideal for ultrasonic cleaning, excellent material compatibility
- Anticorrosive effective



Item	Packaging type	Quantity
2802476	Canister	5L

High Level Instrument Disinfection

- High-level disinfection for all invasive surgical instruments, dental instruments, laparoscopic tools and endoscopes
- Anticorrosive formulation
- Broad-spectrum efficacy: bactericidal, fungicidal, virucidal, and sporicidal



Item	Packaging type	Quantity
2802477	Canister	5L

Instrument Cleaning & Maintenance

Ultrasonic Cleaner

The unit features an automatic locking lid, larger drain, an unique new basket design (hands free basket draining) and touch pad controls.

It comes with an heating and it also features a new easy to clean housing.



Item	Description
2801559	Ultrasonic Cleaner Sweep 2000 with Heat UK Plug
2801560	Ultrasonic Cleaner Sweep 2000 with Heat EU Plug

Basket full size for the Sweep 2000.

Item	Description
2801564	Sweep 2000 Basket Solid with Perforated Bottom 25 x 13 x 7.3 cm

Instrument Disinfection Bath

Instrument bath for disinfecting instruments, including a plastic basket, sieve insert, small lid, and closing lid.

- 4-litre capacity
- Weight: approximately 1.1 kg



Item	Dimensions
2801836	31.5 x 20.6 x 12.5 cm

Complementary Products & Other Accessories

Bin Liner Clear



Item	Dimensions	Quantity
2800452	30 L, 60 x 55 cm, 16 mu	20x20 pcs

Bin Liner Black



Item	Dimensions	Quantity
2800453	30 L, 60 x 55 cm, 31mu	20x20 pcs
2800454	50 L, 75 x 50 cm, 31mu	20x20 pcs
2800455	120 L, 90x65cm, 31mu	20x20 pcs

Room Odour Remover

The fragrant spray creates a pleasant, fresh room fragrance and hides unpleasant odours.



Item	Quantity
2801838	200 ml

Lint Roller Set



Item	Quantity
2800738	1 holder + 2 rolls x 60 sheets

Lint Roller Refill



Item	Quantity
2800739	4 rolls x 60 sheets

Pet Hair Remover

Pet hair remover for dog & cat hair. Fluff & fuzz for the home and car, carpets and upholstery. Lifts, gathers and removes pet hairs.



Item	Quantity
2801010	

covetrus

Care Without Compromise

Products vets trust, value they count on

Hygiene Paper

Paper Towel Roll Midi

Coreless centrefeed towel roll.



Item	Dimensions	Quantity
2800430	2ply, 20 cm x 160 m	6 pcs
2800432	1ply, 20 cm x 280 m	6 pcs

Dispenser

Item	Dimensions
2800422	Midi, White

Paper Towel Roll Mini

Coreless centrefeed towel roll.




Item	Dimensions	Quantity
2802289	2ply, 22 cm x 60 m	12 pcs
2800433	1ply, 20 cm x 120 m	12 pcs

Dispenser

Item	Dimensions
2800423	Mini, White

C-Fold Paper Towels




Item	Dimensions	Quantity
2800424	2ply, 25 x 41 cm, Recycled Paper	20 x 96 pcs
2800426	1ply, 25 x 31 cm,	16 x 152 pcs

Dispenser

Item	Dimensions
2800421	C-fold, Narrow

Z-Fold Towels



Item	Dimensions	Quantity
2800427	1ply, 25 x 23 cm, Recycled Paper	20 x 250 pcs
2800428	2ply, 21 x 24 cm, Z-fold Towels	20 x 160 pcs
2800429	3ply, 32 x 22 cm, Z-fold Towels	20 x 125 pcs

Dispenser

Item	Dimensions
2800420	Z-fold, Wide

Toilet Paper



Item	Dimensions	Quantity
2800434	2ply, 10 cm x 24 m, 48 pcs/ 200 Sheets	bag
2800435	3ply, 10 cm x 31 m, 56 pcs/ 250 Sheets	bag

Facial Tissues Ultrasoft



Item	Dimensions	Quantity
2800830	2ply, 20 x 19.5 cm	100 pcs

Couch Roll Perforated, White

Disposable perforated paper for the exam table. White colour, 2 x 18 grams per m2.



Item	Dimensions	Quantity
2800419	39 cm x 50 m	9 pcs



Personal Protective Equipment (PPE)

Clothes

Disposable Apron PE, Non-sterile



Item	Size	Quantity
2802037	75 x 125 cm	100 pcs
2802038	75 x 140 cm	100 pcs
2802039	75 x 160 cm	100 pcs

Protective Gown with Cuffs PE, Non-sterile



Item	Size	Quantity
2800920	140 cm	1 pc
2800921	156 cm	1 pc
2800922	165 cm	1 pc

Disposable Surgical Gown, Polyethylene

Sterile disposable gowns made of green polyethylene (PE). With tape closures for arms and back. Single use, individually packed.



Item	Colour	Size	Quantity
2800801	●	130 cm	25 pcs

Disposable Overall, Non-woven, Non-sterile



Item	Weight	Size	Quantity
2802908	40g/m ²	M	1 pc
2802909	40g/m ²	L	1 pc
2802910	40g/m ²	XL	1 pc
2802911	40g/m ²	XXL	1 pc

Item	Weight	Size	Quantity
2802912	45g/m ²	M	1 pc
2801483	45g/m ²	L	1 pc
2801484	45g/m ²	XL	1 pc
2801485	45g/m ²	XXL	1 pc

Disposable Surgical Gown, Non-woven, Sterile



Item	Weight	Size	Quantity
2800081	43g/m ²	125 x 150cm	1 pc
2800082	43g/m ²	135 x 155cm	1 pc
2800083	43g/m ²	145 x 164cm	1 pc
2800084	43g/m ²	155 x 164cm	1 pc

Item	Weight	Size	Quantity
2802454	45g/m ²	120 x 140cm	1 pc
2802455	45g/m ²	130 x 150cm	1 pc
2802456	45g/m ²	130 x 160cm	1 pc
2802457	45g/m ²	144 x 170cm	1 pc

Bouffant Cap, Non-woven



Item	Colour	Size	Quantity
2802012	●	61 cm	100 pcs

Shoe Covering

Shoe Covers, Non-woven



Item	Description	Quantity
1	2802013 ● PP, Non-skid, Blue, One size	50 pairs
2	2802027 ● PE, Blue, One size	50 pairs

Disposable Boots

Very strong, non-slip. No welding seams under the foot, comfortable fitting. 100% polyethylen



Item	Quantity
2802916	25 pairs

X-ray Protective Clothing

X-ray Apron Velcro

- With Velcro
- Without pockets
- 0.5 Pb



Item	Description	Size
2802580	● 60 x 100 cm, 0.5 Pb	S
2800827	● 65 x 120 cm, 0.5 Pb	M
2800828	● 75 x 120 cm, 0.5 Pb	L

Thyroid Collar Velcro



Item	Description
2800829	● One Size Blue 0.5 Pb

Latex Examination Gloves

Latex Examination Gloves, Lightly Powdered

- Weight (Grammage/pc): 5.7-6.4
- Protein content (µg/g) ≤ 50
- Thickness Single wall ≥ 0.08 mm

Item	Colour	Size	Quantity
2800109	○	XS	100 / box
2800110	○	S	100 / box
2800111	○	M	100 / box
2800112	○	L	100 / box



Latex Examination Gloves, Powder-free with Aloe Vera

- Weight (Grammage/pc): 5.6-6.2 (size M: 5.8)
- Protein content (µg/g) ≤ 50
- Thickness Single wall ≥ 0.08 mm



Item	Colour	Size	Quantity
2802330	○	XS	100 / box
2802331	○	S	100 / box
2802332	○	M	100 / box
2802333	○	L	100 / box
2802334	○	XL	100 / box

Latex Examination Gloves, Powder-free

- Weight (Grammage/pc): 5.6-6.2 (size M: 5.8)
- Protein content (µg/g) ≤ 50
- Thickness Single wall ≥ 0.08 mm



Item	Colour	Size	Quantity
2802325	○	XS	100 / box
2802326	○	S	100 / box
2802327	○	M	100 / box
2802328	○	L	100 / box
2802329	○	XL	100 / box

Nitrile Examination Gloves

Nitrile Examination Gloves, Powder-free

- Weight (Grammage/pc):
 - All colours: 3.1-3.8
 - Black colour: 3.5-4.3
- Protein content (µg/g) ≤ 50
- Thickness Single wall ≥ 0.05 mm



Item	Colour	Size	Quantity	Item	Colour	Size	Quantity	Item	Colour	Size	Quantity
2800119	○	XS	100 / box	2802335	●	XS	100 / box	2800740	●	XS	200 / box
2800120	○	S	100 / box	2802336	●	S	100 / box	2800741	●	S	200 / box
2800121	○	M	100 / box	2802337	●	M	100 / box	2800742	●	M	200 / box
2800122	○	L	100 / box	2802338	●	L	100 / box	2800743	●	L	200 / box
2800123	○	XL	100 / box	2802339	●	XL	100 / box	2800744	●	XL	200 / box



Nitrile Examination Gloves, Powder-free with Aloe Vera

- Weight (Grammage/pc): 3.1-3.8 (size M: 3.5)
- Protein content (µg/g) ≤ 50
- Thickness Single wall ≥ 0.05 mm

Item	Colour	Size	Quantity
2800129	●	XS	100 / box
2800130	●	S	100 / box
2800131	●	M	100 / box
2800132	●	L	100 / box



Nitrile Examination Gloves, Powder-free, Heavy Duty

- Weight (Grammage/pc): size M: 4.8
- Protein content (µg/g) ≤ 50
- Thickness Single wall ≥ 0.05 mm

Item	Colour	Size	Quantity
2800159	●	XS	100 / box
2800160	●	S	100 / box
2800161	●	M	100 / box
2800162	●	L	100 / box
2800163	●	XL	100 / box



Vinyl Gloves

Vinyl Examination Gloves, Powder-free

- Weight (Grammage/pc): 5.0



Item	Size	Quantity
2800150	S	100 / box
2800151	M	100 / box
2800152	L	100 / box
2800153	XL	100 / box

Gloves Storage

Transparent glove dispenser made of acrylic with recess on the product base for easier removal of the glove boxes. With pre-drilled holes for wall mounting (screws are not included).



Item	Size
2800873	Double
2800874	Single
2800875	Triple

Veterinary Gloves

Veterinary Gloves Supersensitive, Transparent

Fixation band: with elastic tie string for fixation above the elbow: Length: 87 cm, Weight (Grammage/pc): 10 g, Thickness: 0.028 mm


Shoulder protection: the over the head shoulder strap keeps glove in place: Length: 128 cm, Weight (Grammage/pc): 16 g, Thickness: 0.03 mm



Item	Description	Size	Quantity
2800917	With Fixation Band	L	100 / box
2800918	With Fixation Band	S	100 / box
2800919	With Shoulder Protection	One Size	100 / box

Veterinary Gloves Soft

Veterinary gloves made of special soft plastic. One size fits all. No irritation during examination.



Item	Colour	Size	Quantity
2801891	●	Without Shoulder Protection	100 / box
2802035	●	With Shoulder Protection	50 / box

Veterinary Gloves Sensitive

Made of special soft plastic to ensure a supersensitive touch.



Item	Colour	Size	Quantity
2801889	●	Male size	50 pairs
2801890	●	Female size	50 pairs

Surgical Gloves

Latex Surgical Gloves, Powder-free


- Sterile
- Inside surface is coated with an advanced polymer which allow the gloves to be donned easily
- Supersoft latex and anatomical design provide comfort and reduce fatigue
- Microtextured surface gives superior tactile sensitivity in wet and dry conditions



Item	Size	Quantity
2800902	6	50 pairs
2800903	6.5	50 pairs
2800904	7	50 pairs
2800905	7.5	50 pairs
2800906	8	50 pairs
2800907	8.5	50 pairs

Latex Surgical Gloves, Lightly Powdered

- Sterile
- Provide optimal balance between barrier protection and tactile sensitivity
- Supersoft latex film and anatomical design provide unprecedented comfort to reduce fatigue, even during the longest procedures



Item	Size	Quantity
2800910	6	50 pairs
2800911	6.5	50 pairs
2800912	7	50 pairs
2800913	7.5	50 pairs
2800914	8	50 pairs
2800915	8.5	50 pairs

Face Masks

Earloop Mask

- Latex free
- Ear loops are made from special stretch yarn
- 3-ply



Item	Quantity
2802559	50 / box

Particulate Respirator FFP2 Mask



Item	Quantity
2802560	10 / box



Sutures

Suture Overview

Absorbable suture


Product name	Suture material	Suture type	Physical feature	Colour	Coated	Braided	Breaking strength	Absorption rate
PGA	Polyglycolic acid	Synthetic absorbable	Multifilament	Violet	x	x	75% @ 14 days 40% @ 21 days	60-90 days
PGA plus	Polyglycolic acid antibacterial	Synthetic absorbable	Multifilament	Violet	x	x	75% @ 14 days 40% @ 21 days	60-90 days
PGP	PGCL 25 (Polyglecaprone)	Synthetic absorbable	Monofilament	Violet			65% @ 7 days 40% @ 14 days	90-110 days
PDO	Polydioxanone	Synthetic absorbable	Monofilament	Violet			75% @ 14 days 65% @ 28 days 50% @ 42 days	180-220 days
Chromic catgut	Chromic catgut	Natural absorbable	Twisted	Brown			Variable	Generally 21-28 days

Non-absorbable suture

Product name	Suture material	Suture type	Physical feature	Colour	Coated	Braided	Breaking strength	Absorption rate
Polyester green	Polyethylene terephthalate	Synthetic non-absorbable	Multifilament	Green	x	x	N/A	N/A
Polyester white	Polyethylene terephthalate	Synthetic non-absorbable	Multifilament	White	x	x	N/A	N/A
Nylon blue	Dyed Nylon	Synthetic non-absorbable	Monofilament	Blue			N/A	N/A
Copramid	Undyed Nylon	Synthetic non-absorbable	Pseudo-monofilament	White		x	N/A	N/A
Silk	Silk	Natural non-absorbable	Multifilament	Black	x	x	N/A	N/A



Multifilament, braided



Monofilament




Pseudo-monofilament







Twisted multifilament with a monofilament appearance




PGA Braided, coated, polyglycolic acid (PGA), synthetic absorbable suture






Material	Colour(s)	(BSR) Profile	Absorption Rate
PGA	Braided ● Violet	75% @ 14 days 40% @ 21 days	60-90 days

Item	Shape of needle	Needle-code	Needle length mm	EP (metric)	USP	Length in cm	Suture colour	Content per SKU	
2802826	3/8 circle, reverse cutting 	3/8RC19	19	2	3-0	75	Violet	12	
2802827		3/8RC24	24	3	2-0	75	Violet	12	
2802828		3/8RC26	26	2	3-0	70	Violet	12	
2802829		3/8RC30	30	3	2-0	75	Violet	12	
2802830		3/8RC30	30	3.5	0	75	Violet	12	
2802831		3/8RC36	36	4	1	75	Violet	12	
2802816		1/2 circle, reverse cutting 	1/2RC26	26	2	3-0	75	Violet	12
2802817	1/2RC26		26	3	2-0	75	Violet	12	
2802818	1/2RC30		30	3	2-0	75	Violet	12	
2802819	1/2RC30		30	4	1	75	Violet	12	
2802820	1/2RC36		36	3	2-0	90	Violet	12	
2802821	1/2RC36		36	3.5	0	90	Violet	12	
2802822	1/2RC37		37	4	1	90	Violet	12	
2802823	1/2RC40		40	4	1	90	Violet	12	
2802824	1/2RC40		40	5	2	90	Violet	12	
2802825	1/2RC48		48	5	2	90	Violet	12	
2802832	1/2 circle, round bodied, taperpoint 		1/2TP18	18	2	3-0	70	Violet	12
2802833			1/2TP22	22	3	2-0	75	Violet	12
2802834			1/2TP26	26	2	3-0	75	Violet	12
2802835	3/8 circle, round bodied 		1/2TP26	26	3	2-0	75	Violet	12
2802836		1/2TP30	30	3	2-0	75	Violet	12	
2802837		3/8TP24	24	2	3-0	75	Violet	12	
2802838		3/8TP26	26	3.5	0	75	Violet	12	
2802839	without needle	3/8TP32	32	3	2-0	75	Violet	12	
2802840		3/8TP36	36	4	1	75	Violet	12	
2802841					2	3-0	150	Violet	12
2802842					3	2-0	150	Violet	12
2802843					3.5	0	150	Violet	12
2802844					4	1	150	Violet	12
2802845					5	2	150	Violet	12
2802846					6	3+4	150	Violet	12
2802847				7	5	150	Violet	12	


PGA plus Antibacterial-coated, braided polyglycolic acid (PGA), synthetic absorbable suture





Material	Colour(s)	(BSR) Profile	Absorption rate
PGA plus	Braided ● Violet	75% @ 14 days 40% @ 21 days	60-90 days


Item	Shape of needle	Needle-code	Needle length mm	EP (metric)	USP	Length in cm	Suture colour	Content per SKU
2802857	3/8 circle, reverse cutting 	3/8RC19	19	1.5	4-0	70	Violet	12
2802858		3/8RC19	19	2	3-0	45	Violet	12
2802859		3/8RC24	24	2	3-0	70	Violet	12
2802860		3/8RC24	24	3	2-0	70	Violet	12
2802861		3/8RC30	30	3	2-0	70	Violet	12
2802862		3/8RC30	30	3.5	0	70	Violet	12
2802856		1/2 circle, reverse cutting 	1/2RC36	36	4	1	70	Violet
2802863	1/2 circle, round bodied, taperpoint 	1/2TP17	17	1.5	4-0	70	Violet	12
2802864		1/2TP22	22	1.5	4-0	70	Violet	12
2802865		1/2TP22	22	2	3-0	70	Violet	12
2802866		1/2TP22	22	3	2-0	70	Violet	12
2802867		1/2TP26	26	3	2-0	70	Violet	12
2802868	without needle			2	3-0	2x70	Violet	12
2802869				3	2-0	2x70	Violet	12
2802870				3.5	0	2x70	Violet	12
2802871				4	1	2x70	Violet	12
2802872				5	2	2x70	Violet	12



PGP Monofilament poly(glycolide-co-caprolactone) 25 (PGCL) synthetic absorbable suture




Material	Colour(s)	(BSR) Profile	Absorption rate						
PGP (PGCL)	Monofilament ● Violet	65% @ 7 days 40% @ 14 days	90-110 days						
Item	Shape of needle	Needle-code	Needle length mm	EP (metric)	USP	Length in cm	Suture colour	Content per SKU	
2802873	3/8 circle, reverse cutting 	3/8RC19	19	1	5-0	70	Violet	12	
2802874		3/8RC19	19	1.5	4-0	70	Violet	12	
2802875		3/8RC19	19	2	3-0	70	Violet	12	
2802876		3/8RC24	24	2	3-0	70	Violet	12	
2802877		3/8RC24	24	3	2-0	70	Violet	12	
2802878		3/8RC26	26	2	3-0	70	Violet	12	
2802879		3/8RC30	30	3	2-0	70	Violet	12	
2802880		1/2 circle, round bodied, taperpoint 	1/2TP17	17	1.5	4-0	70	Violet	12
2802881			1/2TP22	22	1.5	4-0	70	Violet	12
2802882	1/2TP22		22	2	3-0	70	Violet	12	
2802883	1/2TP26		26	3	2-0	70	Violet	12	



PDO Monofilament polydioxanone (PDO) synthetic absorbable suture




Material	Colour(s)	(BSR) Profile	Absorption rate					
PDO	Monofilament ● Violet	75% @ 14 days 65% @ 28 days 50% @ 42 days	180-220 days					
Item	Shape of needle	Needle-code	Needle length mm	EP (metric)	USP	Length in cm	Suture colour	Content per SKU
2802896	3/8 circle, reverse cutting 	3/8RC12	12	1	5-0	75	Violet	12
2802897		3/8RC19	19	1.5	4-0	75	Violet	12
2802898		3/8RC24	24	2	3-0	75	Violet	12
2802899		3/8RC24	24	3	2-0	75	Violet	12
2802900		3/8RC24	24	3.5	0	75	Violet	12
2802901		1/2 circle, round bodied, taperpoint 	1/2TP22	22	1.5	4-0	75	Violet
2802902	1/2TP26		26	2	3-0	75	Violet	12
2802903	1/2TP30		30	3	2-0	75	Violet	12
2802904	1/2TP36		36	3.5	0	75	Violet	12
2802905	1/2TP48		48	4	1	75	Violet	12



Polyester Braided polyester non-absorbable suture



Material	Colour(s)								
Polyester	Braided ● Green								
Item	Shape of needle	Needle-code	Needle length mm	EP (metric)	USP	Length in cm	Suture colour	Content per SKU	
2802848	3/8 circle, reverse cutting 	3/8RC19	19	1.5	4-0	75	Green	12	
2802849		3/8RC19	19	2	3-0	90	Green	12	
2802850		3/8RC24	24	2	3-0	75	Green	12	
2802851		3/8RC24	24	3	2-0	75	Green	12	
2802852		3/8RC30	30	2	3-0	90	Green	12	
2802853		3/8RC30	30	3	2-0	90	Green	12	
2802854		3/8RC30	30	3.5	0	75	Green	12	
2802855		Straight cutting 	RC48	48	3	2-0	90	Green	12

Nylon/Copramid Nylon monofilament non-absorbable suture




Material	Colour(s)								
Nylon/Copramid	Monofilament/Braided ●○ Blue/White								
Item	Shape of needle	Needle-code	Needle length mm	EP (metric)	USP	Length in cm	Suture colour	Content per SKU	
NYLON									
2802884	3/8 circle, reverse cutting 	3/8RC16	16	1	5-0	75	Blue	12	
2802885		3/8RC19	19	1.5	4-0	75	Blue	12	
2802886		3/8RC19	19	2	3-0	75	Blue	12	
2802887		3/8RC24	24	2	3-0	75	Blue	12	
2802888		3/8RC24	24	3	2-0	75	Blue	12	
2802889		3/8RC30	30	3	2-0	75	Blue	12	
2802890		3/8RC30	30	3.5	0	75	Blue	12	
COPRAMID									
2802893		1/2 circle, round bodied, taperpoint 	1/2TP19	19	2	3-0	75	White	12
2802894	1/2TP24		24	2	3-0	75	White	12	
2802895	1/2TP24		24	3.5	0	75	White	12	
NYLON									
2802891	1/2TP20		20	2	3-0	75	Blue	12	
2802892	1/2TP20	20	3	2-0	75	Blue	12		

Spool Sutures

Absorbable Sutures

Chromic catgut

Twisted gut natural absorbable suture




	Material	Colour(s)	(BSR) Profile	Absorption rate
Chromic catgut	Twisted	● Beige	Variable	Variable (generally 21-28 days)

Item	USP	EP (metric)	Length in m	Colour
2802756	3-0	3	50	Natural brown
2802757	2-0	3.5	50	Natural brown
2802758	0	4	50	Natural brown
2802759	1	5	50	Natural brown
2802760	2	6	45	Natural brown
2802761	3	7	40	Natural brown
2802762	4	8	25	Natural brown



PDO

Monofilament polydioxanone (PDO) synthetic absorbable suture




	Material	Colour(s)	(BSR) Profile	Absorption rate
PDO	Monofilament	● Violet	75% @ 14 days 65% @ 28 days 50% @ 42 days	180-220 days

Item	USP	EP (metric)	Length in m	Colour
2802772	4-0	1.5	25	violet
2802773	3-0	2	25	violet
2802774	2-0	3	25	violet
2802775	0	3.5	25	violet
2802776	1	4	25	violet
2802777	2	5	25	violet
2802778	3+4	6	25	violet



PGA

Braided, coated, polyglycolic acid (PGA), synthetic absorbable suture



	Material	Colour(s)	(BSR) Profile	Absorption rate
PGA	Braided	● Violet	75% @ 14 days 40% @ 21 days	60-90 days


Item	USP	EP (metric)	Length in m	Colour
2802779	4-0	1.5	25	violet
2802780	3-0	2	25	violet
2802781	2-0	3	25	violet
2802782	0	3.5	25	violet
2802783	1	4	25	violet
2802784	2	5	25	violet
2802785	3+4	6	25	violet
2802786	5	7	25	violet
2802787	6	8	25	violet



Non-absorbable Sutures

Copramid

Braided nylon monofilament non-absorbable suture




	Material	Colour
Copramid	Braided	○ White

Item	USP	EP (metric)	Length in m	Colour
2802763	4-0	1.5	75	white
2802764	3-0	2	75	white
2802765	2-0	3	75	white
2802766	0	3.5	75	white
2802767	1	4	75	white
2802768	2	5	50	white
2802769	3+4	6	50	white
2802770	5	7	50	white
2802771	6	8	25	white



Silk

Braided natural silk non-absorbable suture




Material	Colour(s)
Silk Braided	●○ Black/White

Item	USP	EP (metric)	Length in m	Colour
2802788	4-0	1.5	100	black
2802789	3-0	2	100	black
2802790	2-0	3	100	black
2802791	0	3.5	75	black
2802792	1	4	75	black
2802793	2	5	50	black
2802794	3+4	6	50	black
2802795	5	7	50	black
2802796	6	8	25	white
2802797	7	9	25	white
2802798	8	10	25	white



Polyester

Braided polyester non-absorbable suture




Material	Colour(s)
Polyester Braided	○● White/Green

Item	USP	EP (metric)	Length in m	Colour
2802799	2-0	3	100	white
2802800	0	3.5	75	white
2802801	1	4	75	white
2802802	2	5	50	white
2802803	3+4	6	50	white
2802804	5	7	50	white
2802805	6	8	50	white
2802806	7	9	40	white
2802807	8	10	30	white
2802808	4-0	1.5	100	green
2802809	3-0	2	100	green
2802810	2-0	3	100	green
2802811	0	3.5	75	green
2802812	1	4	75	green
2802813	2	5	50	green
2802814	3+4	6	50	green
2802815	5	7	40	green




Suture Needles


Type of circle




Type of eye




Straight



1/2 Circle




3/8 Circle



Covetrus needles are crafted from special hardened, hard-fibre stainless steel to ensure superior sharpness, high ductility, and exceptional bending strength. Precision grinding and dimensional stability guarantee optimal positioning in the needle holder, contributing to a high level of safety and performance in clinical use.




Spring Eye, 12/pk





Item	Shape of Needle	Needle Code	Length mm
2801440	3/8 circle, reverse cutting	3/8RC16	16
2801441		3/8RC20	20
2801442		3/8RC30	30
2801443		3/8RC40	40
2801444		3/8RC45	45
2801445		3/8RC70	70
2801446		3/8RC90	90
2801447	3/8 circle, taperpoint	3/8TP16	16
2801448		3/8TP20	20
2801449		3/8TP30	30
2801450		3/8TP40	40
2801451		3/8TP45	45
2801452		3/8TP70	70
2801453		3/8TP90	90
2801390	1/2 circle, reverse cutting	1/2RC16	16
2801391		1/2RC20	20
2801392		1/2RC30	30
2801393		1/2RC40	40
2801394		1/2RC55	55
2801395	1/2 circle, round bodied, taperpoint	1/2TP16	16
2801396		1/2TP20	20
2801397		1/2TP30	30
2801398		1/2TP40	40
2801399		1/2TP55	55



Regular Eye, 12/Pk










Item	Shape of Needle	Needle Code	Size	Length mm
2801400	3/8 circle, reverse cutting 	3/8RC80	SZ 1	80
2801411		3/8RC75	SZ 2	75
2801413		3/8RC70	SZ 3	70
2801414		3/8RC65	SZ 4	65
2801415		3/8RC60	SZ 5	60
2801416		3/8RC55	SZ 6	55
2801417		3/8RC53	SZ 7	53
2801418		3/8RC50	SZ 8	50
2801419		3/8RC45	SZ 9	45
2801401		3/8RC43	SZ 10	43
2801402		3/8RC41	SZ 11	41
2801403		3/8RC38	SZ 12	38
2801404		3/8RC36	SZ 13	36
2801405		3/8RC33	SZ 14	33
2801406		3/8RC32	SZ 15	32
2801407		3/8RC31	SZ 16	31
2801408		3/8RC28	SZ 17	28
2801409		3/8RC26	SZ 18	26
2801410		3/8RC25	SZ 19	25
2801412		3/8RC24	SZ 20	24
2801420	3/8 circle, taperpoint 	3/8TP80	SZ 1	80
2801431		3/8TP75	SZ 2	75
2801433		3/8TP70	SZ 3	70
2801434		3/8TP65	SZ 4	65
2801435		3/8TP60	SZ 5	60
2801436		3/8TP55	SZ 6	55
2801437		3/8TP53	SZ 7	53
2801438		3/8TP50	SZ 8	50
2801439		3/8TP45	SZ 9	45
2801421		3/8TP43	SZ 10	43
2801422		3/8TP41	SZ 11	41
2801423		3/8TP38	SZ 12	38
2801424		3/8TP36	SZ 13	36
2801425		3/8TP33	SZ 14	33
2801426	3/8TP32	SZ 15	32	
2801427	3/8TP31	SZ 16	31	
2801428	3/8TP28	SZ 17	28	
2801429	3/8TP26	SZ 18	26	
2801430	3/8TP25	SZ 19	25	
2801432	3/8TP24	SZ 20	24	
2801389	1/2 circle, round bodied, taperpoint 	1/2TP31	SZ 16	32

Item	Shape of Needle	Needle Code	Size	Length mm
2801474	Regular eye, 1/2 circle, taperpoint 	1/2TP31	SZ 16	31
2801369	1/2 circle, reverse cutting 	1/2RC80	SZ 1	80
2801380		1/2RC75	SZ 2	75
2801382		1/2RC70	SZ 3	70
2801383		1/2RC65	SZ 4	65
2801384		1/2RC60	SZ 5	60
2801385		1/2RC55	SZ 6	55
2801386		1/2RC53	SZ 7	53
2801387		1/2RC50	SZ 8	50
2801388		1/2RC45	SZ 9	45
2801370		1/2RC43	SZ 10	43
2801371		1/2RC41	SZ 11	41
2801372		1/2RC38	SZ 12	38
2801373		1/2RC36	SZ 13	36
2801374		1/2RC33	SZ 14	33
2801375	1/2RC32	SZ 15	32	
2801376	1/2RC31	SZ 16	31	
2801377	1/2RC28	SZ 17	28	
2801378	1/2RC26	SZ 18	26	
2801379	1/2RC25	SZ 19	25	
2801381	1/2RC24	SZ 20	24	
2801454	Half-Curved Cutting 	1/2C80	SZ 1	80
2801455		1/2C43	SZ 10	43
2801456		1/2C41	SZ 11	41
2801457		1/2C38	SZ 12	38
2801458		1/2C36	SZ 13	36
2801459		1/2C33	SZ 14	33
2801460		1/2C32	SZ 15	32
2801461		1/2C31	SZ 16	31
2801462		1/2C28	SZ 17	28
2801463		1/2C26	SZ 18	26
2801464		1/2C25	SZ 19	25
2801465		1/2C75	SZ 2	75
2801466		1/2C24	SZ 20	24
2801467		1/2C70	SZ 3	70
2801468	1/2C65	SZ 4	65	
2801469	1/2C60	SZ 5	60	
2801470	1/2C55	SZ 6	55	
2801471	1/2C53	SZ 7	53	
2801472	1/2C50	SZ 8	50	
2801473	1/2C45	SZ 9	45	
2801474	Half-Curved Taperpoint 	1/2C31	SZ 16	31

Sutures

Sutures

Item	Shape of Needle	Needle Code	Size	Length mm
2801480		StraightTP38	SZ 12	38
2801481		StraightTP32	SZ 14	32
2801475		StraightRC41	SZ 10	41
2801476		StraightRC38	SZ 12	38
2801477		StraightRC33	SZ 14	33
2801478		StraightRC31	SZ 16	31
2801479		StraightTP41	SZ 10	41

Surgical Skin Adhesive

Surgical Skin Glue

- Fast and painless closure of skin injuries to protect against infections
- Contains n-butyl cyanoacrylate
- Accurate application of the adhesive
- Stops weaker bleeding and supports the healing process
- Acts as a protective barrier against the ingress of foreign substances

Item	Size
2802906	2 ml



Skin Staplers

Skin Staplers, Single Use

Can be used in various surgical procedures for routine skin closure.

- Ensures a perfect wound closure
- Staplers are sterile and ready for use



Item	Description	Size
2802446	35 wide staplers	7.2 mm x 4.9 mm
2802447	35 regular staplers	5.9 mm x 3.9 mm

Single Use Staple Remover

Designed for quick removal of staples without causing damage.



Item	Description
2800577	Staple Remover



Wound Care, Skin Care & Bandages

Cohesive Bandage

Regular Cohesive Bandage

Self-adhering, well flexible and breathable. No scissors needed.



Item	Colour	Dimension	Quantity
2800055	●	5 x 4.5 m	1 roll
2800056	●	5 x 4.5 m	1 roll
2800057	●	5 x 4.5 m	1 roll
2800058	●	5 x 4.5 m	1 roll
2800059	●	5 x 4.5 m	1 roll
2800060	○	5 x 4.5 m	1 roll
2800061	●	5 x 4.5 m	1 roll
2800063	●	7.5 x 4.5 m	1 roll
2800064	●	7.5 x 4.5 m	1 roll
2800065	●	7.5 x 4.5 m	1 roll
2800066	●	7.5 x 4.5 m	1 roll
2800067	●	7.5 x 4.5 m	1 roll
2800068	○	7.5 x 4.5 m	1 roll
2800069	●	7.5 x 4.5 m	1 roll
2800071	●	10 x 4.5 m	1 roll
2800072	●	10 x 4.5 m	1 roll
2800073	●	10 x 4.5 m	1 roll
2800074	●	10 x 4.5 m	1 roll
2800075	●	10 x 4.5 m	1 roll
2800076	○	10 x 4.5 m	1 roll
2800077	●	10 x 4.5 m	1 roll

Vetflex Cohesive Bandage



Item	Colour	Dimension	Quantity
2800185	●	5 x 4.5 m	1 roll
2800186	●	5 x 4.5 m	1 roll
2800187	●	5 x 4.5 m	1 roll
2800188	●	5 x 4.5 m	1 roll
2800189	●	5 x 4.5 m	1 roll
2800191	●	7.5 x 4.5 m	1 roll
2800192	●	7.5 x 4.5 m	1 roll
2800193	●	7.5 x 4.5 m	1 roll
2800194	●	7.5 x 4.5 m	1 roll
2800195	●	7.5 x 4.5 m	1 roll
2800197	●	10 x 4.5 m	1 roll
2800198	●	10 x 4.5 m	1 roll
2800199	●	10 x 4.5 m	1 roll
2800200	●	10 x 4.5 m	1 roll
2800201	●	10 x 4.5 m	1 roll

Cohesive Bandages, Mixed pack with Paw Print

Self-adhering, well flexible and breathable. No scissors needed.

Item	Dimension	Quantity
2800062	5 x 4.5 m (stretched)	36 rolls
2800070	7.5 x 4.5 m (stretched)	24 rolls
2800078	10 x 4.5 m (stretched)	18 rolls
2802430	15 x 4.5 m (stretched)	12 rolls



Care Without Compromise

Products vets trust, value they count on



Orthopaedic Casting

Orthopaedic Casting Tape

Durable fibreglass casting bandages. Stronger than polyester, easy to cut, and more X-ray translucent than plaster. Lightweight, breathable, less sticky, and highly conformable. Becomes rigid about 3–4 minutes after water immersion and supports weight-bearing after 20 minutes. Compatible with cast shoes and bracing devices.



Item	Colour	Width	Length	Quantity
2800582	○	2.5 cm	1.8 m	1 roll
2800583	○	5.0 cm	3.6 m	1 roll
2800584	○	7.5 cm	3.6 m	1 roll
2800585	○	10 cm	3.6 m	1 roll
2800586	●	2.5 cm	1.8 m	1 roll
2800587	●	5.0 cm	3.6 m	1 roll
2800588	●	7.5 cm	3.6 m	1 roll
2800589	●	10 cm	3.6 m	1 roll
2800590	●	5.0 cm	3.6 m	1 roll
2800591	●	7.5 cm	3.6 m	1 roll
2800592	●	10 cm	3.6 m	1 roll

Bandage Covering

Abdominal Bandage

Protects the surgery wound and offers protection against skin issues. Good acceptance by the animals also for longer term use. Washable frottee material made of 100% cotton and adjustable through Velcro.



Item	Size	Length
2800668	Size 1 for Cats	33 cm
2800669	Size 2 for Poodles	36 cm
2800670	Size 3 for Spaniel	44 cm
2800671	Size 4 for Small German Shepherds	49 cm
2800672	Size 5 for German Shepherds	55 cm
2800673	Size 6 for Great Dane	57 cm
2800674	Size 7 for Saint Bernard's	64 cm

Tubular Bandage

Seamless, breathable, and highly elastic. Made of 67% cotton and 33% rayon for comfortable wear. Ideal for skin protection, under casts, adhesive or ointment dressings, and under compression bandages.



Item	Size	Length	Quantity
2800834	Size 0, (flat) width 1.5 cm	20 m	1 roll
2800835	Size 1, (flat) width 2.5 cm	20 m	1 roll
2800836	Size 2, (flat) width 6.0 cm	20 m	1 roll
2800837	Size 3, (flat) width 8.0 cm	20 m	1 roll
2800838	Size 4, (flat) width 10 cm	20 m	1 roll
2800839	Size 5, (flat) width 12 cm	20 m	1 roll
2800840	Size 6, (flat) width 16 cm	20 m	1 roll
2800841	Size 7, (flat) width 21 cm	20 m	1 roll

Tape

Adhesive Wound Dressing on a Roll

Flesh-colored, non-elastic cotton tape with strong zinc oxide adhesive. Easy to tear, breathable, and writable. Ideal for securing dressings, tubes, cannulae, catheters, and probes.



Item	Color	Dimensions	Quantity
2800842	●	5 m x 1.25 cm	18 rolls
2800843	●	5 m x 2.50cm	12 rolls
2800844	●	5 m x 5 cm	6 rolls

Silk Adhesive Wound Dressing on a Roll

White cellulose acetate with skin-friendly adhesive. Easy to tear, breathable, and writable. Ideal for securing dressings, tubes, cannulae, catheters, and probes.



Item	Quantity	Dimensions	Quantity
2800845	○	5 m x 1.25 cm	18 rolls
2800846	○	5 m x 2.50cm	12 rolls
2800847	○	5 m x 5 cm	6 rolls

Non-woven Paper Tape

This white paper tape is made from non-elastic, micro-perforated fleece with a skin-friendly polyacrylate adhesive. It can be torn both crosswise and lengthwise, is breathable, vapour-permeable, and can be written on with a ballpoint pen. It is suitable for fixing dressings and securing tubes, cannulas, catheters, or probes.



Item	Dimensions	Quantity
2800848	5 m x 1.25 cm	18 rolls
2800849	5 m x 2.50 cm	12 rolls
2800850	5 m x 5 cm	6 rolls

Transparent Adhesive Tape

This transparent tape is made from perforated polyethylene film with a skin-friendly polyacrylate adhesive. It can be torn crosswise, lengthwise, or down the middle, and is breathable, vapour-permeable, and easy to write on with a ballpoint pen. Ideal for fixing dressings and securing tubes, cannulas, catheters, or probes.



Item	Dimensions	Quantity
2800854	5 m x 1.25 cm	18 rolls
2800855	5 m x 2.50cm	12 rolls
2800856	5 m x 5cm	6 rolls

Sport Tape

This non-elastic tape made of 100% cotton is coated with an extra strong zinc oxide caoutchouc. It adheres safely on the skin as well as on under-tape bandages.



Item	Dimensions	Quantity
2800851	10 m x 2.5 cm	1 roll
2800852	10 m x 3.8 cm	1 roll
2800853	10 m x 5 cm	1 roll

Adhesive Fibre Ribbon Tape

Item	Description	Dimensions
2800172	● Black - Tearable	50 m x 1.9 cm
2800173	● Black - Tearable	50 m x 3.8cm
2800175	● Black - Tearable	50 m x 6.0cm
2800176	● Black	50 m x 5.0cm
2800177	○ Off White	50 m x 1.9 cm
2800178	○ Off White	50 m x 2.5 cm
2800179	○ Off White	50 m x 3.8cm
2800180	○ Off White	50 m x 5.0cm
2800181	○ Waterproof, White	50 m x 1.9 cm
2800182	○ Waterproof, White	50 m x 2.5 cm
2800183	○ Waterproof, White	50 m x 3.8cm
2800184	○ Waterproof, White	50 m x 5.0cm



Padding

VET Padding Bandage

Unbleached polyester bandage for padding and stuffing. This VET bandage is thicker and slightly stiffer than the regular 'Synthetic padding bandage' 2800857-2800859.



Item	Dimensions	Quantity
2800513	3 m x 6 cm	1 roll
2800514	3 m x 10 cm	1 roll
2800515	3 m x 15 cm	1 roll
2800516	3 m x 25 cm	1 roll

Synthetic Padding Bandage

This synthetic padding is made from non-woven polyester filaments, mechanically bonded without chemical agents. It is ideal for use under rigid and semi-rigid bandages—such as casts or zinc adhesive bandages—and for padding splints or cushioning compression bandages. Softer and more flexible than the VET variant (2800513 to 2800516).



Item	Dimensions	Quantity
2800857	3 m x 10 cm	72 rolls
2800858	3 m x 15 cm	36 rolls
2800859	3 m x 6 cm	108 rolls

Cotton Wool Roll, 100% Cotton



Item	Dimensions	Thickness	Quantity
2801955	20 cm x 7 m	675 g/m ²	1000 gr
2801956	20 cm x 3.5 m	500 g/m ²	500 gr

Cotton Wool Roll – Premium Quality

Premium-quality, highly absorbent cotton wool roll made of 95% cotton and 5% viscose. The 500 g roll is pH balanced, naturally white without harmful chemicals, and offers reliable absorbency. It features a consistent layer, clean-cut edges, and a very soft texture.



Item	Description	Quantity
2802022	Absorbent 380 g/m ²	500 g

Medical bleached absorbent cotton wool

High-quality cotton wool covered in gauze, ideal for wound protection and management. It insulates and cushions wounds, can be applied over dressings before bandaging, and serves as padding for hoof poultices. Also suitable as a compression pad for bleeding wounds and for general swabbing and cleaning.



Item	Dimensions	Quantity
2802236	9 cm x 6 m	200 g/roll
2802237	18 cm x 6 m	340 g/roll
2802238	30 cm x 3.2 m	500 g/roll
2802239	45 cm x 2.3 m	500 g/roll

Cotton Wool – Zigzag, 100% Cotton

Item	Quantity
2800555	100 g
2800556	250 g
2800557	500 g



Other Bandages

Teat Bandage

A unique, elastic, absorbent, and cohesive bandage that adheres only to itself – not to skin, hair, or wounds. It is water-resistant, latex-free, and compatible with ointments and creams. The bandage is soft, easy to apply, and remains fully flexible.



Item	Dimensions	Quantity
2802449	6 cm x 5 m	1 roll

Elastic Adhesive Bandage

100% cotton bandage, gorge edged and with natural rubber coating.



Item	Dimensions	Quantity
2800051	4.5 m x 2.5 cm (1")	1 roll
2800052	4.5 m x 5.0 cm (2")	1 roll
2800053	4.5 m x 7.5 cm (3")	1 roll
2800054	4.5 m x 10 cm (4")	1 roll

Elastic Conforming Bandage, Single Packed

Woven edge, 45% cotton and 55% polyamide.



Item	Dimensions	Quantity
2800517	4 cm x 4 m	20 rolls
2800518	6 cm x 4 m	20 rolls
2800519	8 cm x 4 m	20 rolls
2800520	10 cm x 4 m	20 rolls
2800521	12 cm x 4 m	20 rolls

Vet Polster Bandages

Elastic, self-adhesive bandage which offers excellent support and comfort. It is thicker (0.5 cm) and more robust than standard teat bandages. It is easy to apply, does not stick to skin, hair, or wounds, and adapts to any body part.



Item	Description	Quantity
2802654	Vet Polster Bandage (Reinforced)	1 roll
2802655	Vet Polster Bandage	1 roll

Cotton Bandage

These non-sterile bandages are primarily used to protect dressings and keep them securely in place. Made from 100% non-bleached cotton, they offer a reliable grip between layers and help secure splints effectively.



Item	Dimensions	Quantity
2800897	4 m x 4 cm	20 rolls
2800898	4 m x 6 cm	20 rolls
2800899	4 m x 8 cm	20 rolls
2800900	4 m x 10 cm	20 rolls

Ideal Bandage

This air-permeable, permanently elastic, and durable bandage offers a secure, long-lasting fit and lies close to the skin. It is ideal for fixing dressings and providing support or relief in primary care. The material consists of 67% cotton, 28% polyamide, and 5% latex.



Item	Dimensions	Quantity
2800878	5 m x 4 cm	10 rolls
2800879	5 m x 6 cm	10 rolls
2800880	5 m x 8 cm	10 rolls
2800881	5 m x 10 cm	10 rolls
2800882	5 m x 12 cm	10 rolls

Non-woven

Non Woven Sponge, Non-sterile

Made from synthetic polyester.



Item	Dimensions	Ply	Quantity
2800214	5 x 5 cm	4ply	200 / pack
2802957	7.5 x 7.5 cm	4ply	200 / pack
2800216	10 x 10 cm	4ply	200 / pack

Dressing Gauze, 17-thread, Sterile

Made from 100% cotton, this individually packed gauze offers high absorption for various procedures, reducing the number of swabs needed and improving cost-efficiency. Folded edges and stable threads prevent fiber residue in wounds, supporting faster healing. Suitable for wound dressing, post-surgery coverage, protection, cleaning, and incision drainage.



Item	Dimensions	Quantity
2802451	0.25 m ²	1pc
2802452	0.50 m ²	1pc
2802453	1.00 m ²	1pc

Gauze Swabs, Non-sterile – 8 Layer

Made from synthetic polyester.



Item	Dimensions	Ply	Quantity
2800229	7.5 x 7.5 cm	8ply	2500 / carton

Hydrophilic

Gauze Swabs Hydrophilic, Non-sterile

100% cotton, edge folded and highly absorbent.



Item	Dimensions	Ply	Quantity
2800217	5 x 5 cm	12ply	100 / pack
2800218	5 x 5 cm	8ply	100 / pack
2800219	7.5 x 7.5 cm	12ply	100 / pack
2800220	7.5 x 7.5 cm	8ply	100 / pack
2800221	10 x 10 cm	12ply	100 / pack
2800222	10 x 10 cm	8ply	100 / pack

Gauze Swabs Hydrophilic, Sterile

Edge folded and highly absorbent.



Item	Dimensions	Ply	Quantity
2802965	5 x 5 cm	12ply	2 pcs / pouch 100 pcs
2800225	7.5 x 7.5 cm	12ply	2 pcs / pouch 80 pcs
2802967	10 x 10 cm	12ply	2 pcs / pouch 50 pcs

Gauze Swabs X-ray Hydrophilic Sterile

Sterile hydrophilic gauze compresses with X-ray contrast thread.



Item	Dimensions	Ply	Quantity
2800226	5 x 5 cm	12ply	10 x 10 pcs
2802969	7.5 x 7.5 cm	12ply	10 x 10 pcs
2800228	10 x 10 cm	12ply	10 x 10 pcs

Gauze & Cotton Balls

Gauze Balls, Non-sterile

100% Cotton, fibre free and non-chlorine bleached.

Item	Dimensions	Size	Quantity
2800507	6 x 9 mm	Hazelnut size	500 / pack
2800508	15 x 15 mm	Walnut size	500 / pack
2800509	20 x 20 mm	Plum size	500 / pack
2800510	25 x 25 mm	Egg size	500 / pack
2800511	30 x 30 mm	L	250 / pack
2800512	30 x 40 mm	XL	250 / pack



Cotton Wool Balls, Non-sterile, 100% Cotton



Item	Quantity
2800558	250 / pack

Cellulose Compress

Compresses are bleached. Ideal for cleaning & wiping of the skin.



Item	Size	Quantity
2800169	4 x 5 cm	2 rolls

Cellulose Paper Wipes

Non-bleached, cutted into 2 sizes.



Item	Size	Quantity
2800170	20 x 30 cm	3 x 5 kg
2800171	40 x 60 cm	3 x 5 kg

Wound Care

Manuka Honey Gel

Sterile medical-grade Manuka honey helps maintain a moist wound environment that supports healing.

Rich in Methylglyoxal (MGO), it is standardized for antibacterial efficacy and sterilized for safety. Sourced from 100% Leptospermum scoparium honey from New Zealand, it promotes autolytic debridement, supports granulation, and reduces wound odor. Its low pH prevents bacterial growth.

Under veterinary supervision, it can be used for trauma, surgical, ulcerated, infected, and persisting wounds, as well as 1st and 2nd degree burns, lacerations, and abrasions. Available in 15 g and 42.5 g single-use tubes, it minimizes cross-contamination and can be applied directly to the wound or dressing.

Item	Quantity
2800417	15 g
2800418	42.5 g



Skin Care

Blue Spray

Antibiotic-free blue spray for skin and claw. It offers protection and care for severely stressed and sensitive areas of skin.



Item	Size
2801498	200 ml

Zincoxide Spray

Ointment spray that supports natural skin regeneration in animals. The formula includes zinc oxide, a key component for skin regeneration.



Item	Size
2801501	200 ml

Powder Spray

Antibiotic-free microfine powder spray designed to protect and care for sensitive or strained skin areas in animals. It supports natural regeneration of skin and fur, and the powder helps dry out the sprayed area.



Item	Size
2801500	200 ml

Aluminium Spray Silver Coloured

The spray creates a film for the protection of sensitive skin against dirt and other environmental influences.



Item	Size
2801499	200 ml

Adhesive Remover Spray

Removes adhesives and therefore makes the removal of bandage and tape faster, easier and more comfortable. The Spray is formulated to reduce stress and discomfort during dressing.



Item	Size
2802656	50ml
2802657	200ml

Miscellaneous

Bitter Paste / Spray

The items have an extremely bitter taste for animals (cats, dogs, horses) and are ideal for preventing licking and biting on surfaces and bandages. It contains natural pepper extract (capsaicin) as an active ingredient.



Item	Description	Size
2801502	Paste	250 ml
2801490	Spray	100 ml

Hot/Cold

Hot/Cold Packs – Reusable

Great for treating a wide range of soft tissue injuries, Pack will remain flexible at low temperatures.



Item	Size	Quantity
2800876	14 x 13 cm	10 pcs
2800877	29 x 12 cm	10 pcs



Pharmaceuticals & Chemicals



Some products from this chapter may be available only in selected countries due to the local registration/regulatory requirements. Please check availability with your Covetrus representative.

Prescription Vials

Prescription Vial with Screw Cap, Child Resistant

Item	Description	Size
2802040	6 dram	24 ml
2802041	8 dram	30 ml
2802042	13 dram	50 ml
2802043	16 dram	60 ml
2802044	20 dram	75 ml
2802045	30 dram	120 ml
2802046	40 dram	160 ml



Pharmaceuticals

Isoflurane - Inhalation Vapour

Inhalation anaesthetic for veterinary use with active substance Isoflurane (1000 mg/g). Isoflurane is always administered in conjunction with air and/or pure oxygen. Isothesia is liquid at room temperature.

- To be used for dogs, cats, ferrets, horses, reptiles, ornamental birds and small mammals
- Compatible with all isoflurane specific vapourisers
- Prescription only medicine

Item	Description	Quantity
2800026	Isothesia - Inhalation Vapour, liquid, Germany	250 ml
2800028	Isothesia - Inhalation Vapour, liquid, Ireland	250 ml
2800029	Isothesia - Inhalation Vapour, liquid, Poland	250 ml
2800030	Isothesia - Inhalation Vapour, liquid, Netherlands	250 ml
2800032	Isothesia - Inhalation Vapour, liquid, Greece	250 ml
2800033	Isothesia - Inhalation Vapour, liquid, Romania	250 ml
2800691	Isothesia - Inhalation Vapour, liquid, Belgium	250 ml
2800024	Isothesia NXT- Inhalation Vapour, liquid, Australia	250 ml
2800692	Isothesia NXT - Inhalation Vapour, liquid, New Zealand	250 ml
2802569	Isothesia - Inhalation Vapour, liquid, Czech Republic	250 ml
2800025	Isofane - Inhalation Vapour, liquid, UK	250 ml
2800027	Isofane - Inhalation Vapour, liquid, Spain	250 ml
2802732	Attane - Inhalation Vapour liquid, Switzerland	250 ml



Tulathromycin - Solution for injection

Fast-acting, broad-spectrum macrolide antibiotic, designed for the treatment and metaphylaxis of respiratory diseases in cattle, pigs, and sheep. It is also indicated for infectious bovine keratoconjunctivitis and early-stage foot rot in sheep.

- Reliable respiratory protection: Targets key pathogens like Mannheimia haemolytica, Pasteurella multocida, and Mycoplasma bovis
- Flexible administration: Subcutaneous use in cattle; intramuscular use in pigs and sheep.
- Prescription only medicine



Ireland only

Item	Description	Quantity
2803083	Latroxin - Solution for Injection, 100mg/ml	100 ml
2803095	Latroxin - Solution for Injection, 100mg/ml	250 ml

Amoxicillin and Clavulanic Acid - Chewable Tablets

Flavoured and chewable tablets combining Amoxicillin and Clavulanic Acid. For use in dogs and cats for the treatment of bacterial infections affecting the respiratory, digestive, or urinary systems, as well as skin and dental infections.

- One side scored for easy and accurate dosing
- Chewable formulation for easier application
- Enhanced flavour and palatability
- Prescription only medicine



Ireland only

Item	Description	Quantity
2802566	CLAMADEN - Chewable Tablets, 40mg/10mg	100 Tabs
2802567	CLAMADEN - Chewable Tablets, 200mg/50mg	100 Tabs

Carprofen - Tablets

Exclusively for canines for the reduction of inflammation and pain caused by musculoskeletal disorders and degenerative joint disease. As a follow up to parenteral analgesia in the management of post operative pain.

- Round, dark brown, marbled tablets with visible darker spot
- One side scored for easy and accurate dosing
- Enhanced flavour and palatability
- Prescription only medicine



Ireland only

Item	Description	Quantity
2803096	KARPROVET - Tablets, 20mg	100 Tabs
2802565	KARPROVET - Tablets, 50mg	100 Tabs
2803097	KARPROVET - Tablets, 100mg	100 Tabs

Carprofen - Solution for Injection

Non-steroidal anti-inflammatory drug (NSAID) designed for the control of post-operative pain and inflammation in dogs and cats.

Dogs: For the control of post-operative pain and inflammation following orthopaedic and soft tissue (including intraocular) surgery.
Cats: For the control of post operative pain following surgery.

- Clear, pale yellow coloured solution
- Prescription only medicine

Ireland only

Item	Description	Quantity
2802564	KARPROVET - Solution for Injection, 50mg/ml	20 ml



Meloxicam range

Comprehensive range of 10 formulations including oral suspensions for cats, guinea pigs, and dogs, tablets for cats, and injectable solutions for cattle, pigs, horses, dogs, and cats – designed to provide effective pain and inflammation relief across species.

Broad range of presentations

- Oral suspensions
- Injectable solutions
- Chewable tablets
- Effective pain and inflammation relief
- Flexible dosing and administration
- Prescription only medicine

Belgium only

Item	Description	Quantity
2802722	Metaxx 0,5 mg/ml, Oral Suspension	5ml
2802723	Metaxx 0,5 mg/ml, Oral Suspension	10ml
2802724	Metaxx 0,5 mg/ml, Oral Suspension	25ml
2802726	Metaxx 1,5 mg/ml, Oral Suspension	25ml
2802727	Metaxx 1,5 mg/ml, Oral Suspension	50ml
2802728	Metaxx 1,5 mg/ml, Oral Suspension	125ml
2802729	Metaxx 5 mg/ml, Solution for injection	20ml
2802730	Metaxx 5 mg/ml, Solution for injection	100ml
2802731	Metaxx 20 mg/ml, Solution for injection	100ml
2802732	Metaxx 0,25 mg, Tablets	9x10 tablets



Pentobarbital - Solution for Injection

Solution for injection indicated for euthanasia is available in concentrations of 40%. The solution for injection is now licensed for a wide range of target animals, making it even more widely applicable in practice.

Broad range of presentations

- Licensed for 20 target animals in total
- Prescription only medicine

Poland only

Item	Description	Quantity
2803100	Euthanimal 400 mg/ml, Solution for Injection	100 ml



Miconazole / Prednisolone / Polymyxin B - Ear Drops

Formulated for the treatment of primary and secondary skin infections and skin adnexa in dogs, cats, and guinea pigs. Also indicated for otitis externa in dogs and cats caused by miconazole- and polymyxin B-susceptible bacteria and fungi.

Ingredients:

- Miconazole nitrate 23.00 mg
- Prednisolone acetate 5.00 mg
- Polymyxin B sulfate 5500 IU
- Treating otitis externa and skin problems
- Prescription only medicine

Ireland and Greece only

Item	Description	Quantity
2802925	Otomicol - Oral Suspension	15 ml



Chemicals - Water

Purified Water

Suitable for many autoclave sterilisation device, instrument cleansing processors and ultrasonic water bath agitation cleaning devices.

Uk, Ireland only



Item	Size
2800602	5 L
2800603	25 L

Distilled & De-mineralised Water

De-mineralised and lab-tested for quality. Suitable for laboratory use



Item	Description	Size
2802269	Double Distilled	1 L
2801365	Distilled	5 L
2801366	Distilled	4 x 5 L
2801367	Distilled	10 L

Formal Saline (formaldehyde)

A ready to use 10% Formal Saline with formaldehyde and sodium chloride.

Uk, Ireland only



Item	Size
2800608	2 L

Propylene Glycol

Veterinary Grade Propylene Glycol (1,2-Propanediol)

UK, Ireland only



Item	Size
2800618	2 L
2800619	5 L
2800620	25 L

Liquid Paraffin

Veterinary grade Liquid Paraffin.

Uk, Ireland only



Item	Size
2800612	2 L

Iodine

Spirit based. Superior quality. For Drying wounds and umbilicals.

Health & Safety: Highly flammable liquid and vapour. Harmful if inhaled or swallowed, may cause damage to organs.

Causes skin irritation, serious eye irritation. Wear protective gloves, clothing and eye protection.

Active ingredient: Ethanol, Iodine, Potassium Iodide, Methanol, Amines, C10-C-16 Alkyl Derivatives reaction products with Chloroacetic Acid.

Storage: Keep away from heat sources, sparks and open flame and oxidising materials. Keep container tightly closed and sealed until ready for use. Store in original container protected from direct sunlight in a dry, cool and well-ventilated area.

Disposal: This material and its container must be disposed of in a safe way. Dispose of contents/containers to an approved waste disposal plant.

Uk, Ireland only



Item	Description	Size
2800625	Iodine Strong 10% with Sprayhead	500 ml
2800626	Iodine Strong 10%	500 ml
2800627	Iodine Strong 10%	2 L



Some products from this chapter may be available only in selected countries due to the local registration/regulatory requirements. Please check availability with your Covetrus representative.

Chemicals - Alcohol

Denatured Alcohol

Denatured Alcohol

Active ingredient: 94% Ethanol, 3% Isopropanol, 3% Methyl Ethyl Ketone

Uk, Ireland only



Item	Size
2800615	2 L

BP Surgical Alcohol 95%

Pre-operative British Pharmacopeia (BP) grade formulation.

Active ingredient: Ethanol, Methanol, Castor Oil, Diethyl Puthalate, Pyroigneous Acids, Methyl Salicylate

Uk, Ireland only



Item	Size
2800616	1 L
2800617	5 L

Hydrogen Peroxide 20V 6%

A ready to use Hydrogen Peroxide. Active ingredient: Hydrogen Peroxide 6%

Uk, Ireland only



Item	Size
2800610	1 L
2800611	2 L

Copper Sulphate

Active ingredient: Copper Sulphate

Uk, Ireland only



Item	Size
2800604	1 kg

Magnesium Sulphate

Active ingredient: Epsom Salts

Uk, Ireland only



Item	Size
2800605	1 kg
2800606	2 kg

Glucose Powder

Active ingredient: D-Glucose

Uk, Ireland only



Item	Size
2800621	1 kg



Urology

Urinary Catheters

Cat Urinary Catheter

Made out of polypropylene, flexible and with rounded end. Ideal for urine sample collection, urine retention relief and bladder lavage.

- Luer fitting
- Atraumatic side holes (closed end)



Item	Description
2800035	without Stylet 3FG (1.00 x 130 mm)
2800036	without Stylet 4FG (1.33 x 130 mm)
2800037	with Stylet 3FG (1.00 x 130 mm)
2800038	with Stylet 4FG (1.33 x 130 mm)

Dog Urinary Catheter

Made out of polyethylene, flexible and with rounded tip. Ideal for urine sample collection, urine retention relief and bladder lavage.

- Female luer lock fitting
- Atraumatic side holes (closed end)



Item	Description
2800039	10FG (3.33 x 500 mm)
2800040	4FG (1.33 x 500 mm)
2800041	6FG (2.00 x 500 mm)
2800042	8FG (2.66 x 500 mm)

Urine Collection Dogs

Uricollect Dogs

Urine collection device to easily collect urine from bitches, consisting of a transparent sample tube that is attached to a scoop-shaped collection container. The urine automatically flows into the collection tube, which can then be closed with a screw cap for transport.



Item	Description
2800802	Uricollect Urine Collection Device Dog

Urine Collection Cats

Urigrit - Polypropylene Grit

Urigrit is a specially developed polypropylene grit to be used to collect urine from cats. Does not interact with the urine.



Item	Description
2800001	UriGrit Urine Collection Cat

Hydrosand Urine Collection Kit

Hydrophobic litter sand for collecting urine from cats. Its structure resembles bentonite litter.



Item	Description	Quantity
2800890	Kit excl. Tray	3 sets
2800891	Kit excl. Tray	1 set

Urine Container

Urine Cup – Screw Type

Transparent sample cup with 20 ml scale marking. Matching green screw cap (Art. 2800895) available separately.



Item	Description	Size
2800893	Urine Cup Screw Type – excl. Lid	100 ml
2800895	Urine Lid Green for Cup Screw Type	100 ml

Urine Cup – Snap on Type

Transparent sample cup with 25 ml scale marking. Matching green screw cap (Art. 2800896) available separately.



Item	Description	Size
2800894	Urine Cup Snap on Type – excl. Lid	125 ml
2800896	Urine Lid Red for Snap on Cup	125 ml

Fluid Collection Vial

Polypropylene vial with lid and 16 ml capacity, 10 cm length. Non-sterile.



Item	Description	Size	Quantity
2800003	With Lid Polypropylene	16 ml	20 pcs

Diagnostics

Urine Test Strips

Urine Test Strips

Reagent test strips for quick detection of different indications in the urine of small animals. Dip and Read Technique.

- Item 2800567 can also be used with the One Step Plus Urine Analyser

Item	Description	Quantity
2800563	Vet Urine Test Strips Glucose	100 strips
2800565	Vet Urine Test Strips 9 + Leuko	100 strips
2800566	Vet Urine Test Strips 10 Way	100 strips
2800567	Vet Urine Test Strips 11 Way	100 strips



Faeces Analysis

Stool Sample Tube

Dimensions: 22 x 63 mm (ØxL)



Item	Description	Quantity
2800892	Stool Sample Tube	100 pcs

Flotation Liquid

Flotation liquid for fecal analysis. Magnesium sulphate solution.



Item	Description	Size
2801491	Fecalizer Flotation Liquid	1L

Fecal Diagnostic Kit

A disposable fecal collection device designed for a clinic or at home setting. Designed for use as both a fecal collection container and a flotation exam container. Various gastrointestinal parasite eggs can be detected. For use in feces of dogs and Cat.



Item	Description	Quantity
2800168	Fecal Diagnostic Flotation Kit	50 / box

Dermatophyte Tests



Petri Dish with DTM

The fungal culture test consists of individually petri dishes with removable lids that allow for easy inoculation and collection of samples for identification. It contains Dermatophyte Test Medium (DTM) and is designed to diagnose the most relevant veterinary pathogenic fungi called dermatophytes: Microsporium, Trichophyton and Epidermophyton spp. DTM contains agents that suppress the growth of saprophytic fungi and contaminating bacteria, while allowing the growth of dermatophytic fungi. These fungi will turn the colour of DTM from orange to red.

Item	Description	Quantity
2800579	Dermatophyte Test Petri Dish DTM	6 plates

Vial with DTM

This fungal culture consists of a bottle with a removable lid and a plate in it on which fungi can grow. Easily applicable way to inoculate fungi and collect sample for diagnostics. The fungal culture test consists of vials with a slant medium inside, of which the cap can be removed allowing inoculation and collection of samples for identification. Easily applicable way to inoculate fungi and collect sample for diagnostics. It contains Dermatophyte Test Medium (DTM) and is designed to diagnose the most relevant veterinary pathogenic fungi called dermatophytes: Microsporium, Trichophyton and Epidermophyton spp. DTM contains agents that suppress the growth of saprophytic fungi and contaminating bacteria, while allowing the growth of dermatophytic fungi. These fungi will turn the colour of DTM from orange to red.



Item	Description	Quantity
2800580	Dermatophyte Test Vial DTM	10 vials

Three-chambered Petri Dish with DTM / ESA / SAB

The fungal culture test consists of individually wrapped three-chambered petri dishes with removable lids that allow for easy inoculation and collection of samples for identification. Each plate contains a Dermatophyte Test Medium (DTM), an Enhanced Sporulation Agar (ESA) and a Sabouraud agar (SAB) and is designed to diagnose the most relevant veterinary pathogenic fungi called dermatophytes: Microsporium, Trichophyton and Epidermophyton spp. It contains Dermatophyte Test Medium (DTM) and is designed to diagnose the most relevant veterinary pathogenic fungi called dermatophytes: Microsporium, Trichophyton and Epidermophyton spp. Like DTM, ESA contains supplements to inhibit growth of bacteria and saprophytic fungi, but the colour change is from yellowish to blue-green. What makes ESA unique is that it enhances both pigmentation and sporulation of dermatophytes, thus allowing for proper identification of the fungal isolates.



Item	Description	Quantity
2800581	Dermatophyte Test Petri Dish DTM / ESA / SAB	4 plates

Rapid Tests

Covetrus offers a complete range of rapid immunochromatic tests (also known as lateral flow assays). Boxes of 2 tests and 10 tests are available

- User-friendly, immediate on-site diagnostics
- Quick results visible within a few minutes
- Can be stored at room temperature
- No need to wait for them to warm up before use



Overview of available Rapid test

2 tests	10 tests	Antigen (Ag) or Antibody (Ab)	Species	Sample	
2800014	2800004	Rapid Test Feline Leukaemia Virus (FeLV)	Ag	Cat	
2800015	2800005	Rapid Test Feline Immunodeficiency Virus (FIV)	Ab	Cat	
2800016	2800006	Rapid Test Feline Combi (FeLV + FIV)	Ag + Ab	Cat	
-	2800013	Rapid Test Canine Leishmania	Ab	Dog	Blood
-	2802015	Rapid Test Dirofilaria Immitis	Ag	Dog	
2802016	2802017	Rapid Test Dirofilaria Ag + Anaplasma + Ehrlichia + Leishmania Ab Combi	Ag / Ab	Dog	
	2800008	Rapid Test Corona Test (CCV + FCoV)	Ag	Cat & Dog	
2800019	2800009	Rapid Test Giardia	Ag	Cat & Dog	
	2802736	Rapid Test Giardia + CPE	Ag	Cat & Dog	
-	2800010	Rapid Test Parvo + Corona Combi	Ag	Cat & Dog	Faeces
2800021	2800011	Rapid Test Parvo + Corona + Giardia Combi	Ag	Cat & Dog	
2800017	2800007	Rapid Test Parvovirus (CPV + FPV)	Ag	Dog	
	2802737	Rapid Test Rotavirus + Bovine Coronavirus + E.coli K99 + Cryptosporidium parvum	Ag	Bovine	
2800022		Rapid Test Lyme (Borrelia)	Ag	Cat & Dog	Ticks

Overview Sensitivity and Specificity

Product	Item code	Pathogen 1		Pathogen 2		Pathogen 3		Pathogen 4	
		Sensitivity*	Specificity*	Sensitivity*	Specificity*	Sensitivity*	Specificity*	Sensitivity*	Specificity*
FeLV Ag	2800004	94,59%	99,99%						
	2800014	94,59%	99,99%						
FIV Ab	2800005	91,67%	97,83%						
	2800015	91,57%	97,83%						
FeLV Ag + FIV Ab	2800006	94,59%	99,99%	91,67%	97,83%				
	2800016	94,59%	99,99%	91,67%	97,83%				
FCoV/CCoV Ag	2800008	99,99%	97,50%						
	2800018	99,99%	97,50%						
	2800007	99,99%	97,50%						
	2800017	99,99%	97,50%						
Giardia	2800009	96,00%	97,22%						
	2800019	96,00%	97,22%						
Giardia + CPE Ag	2802736	96,00%	97,22%	97,83%	98,15%				
	2800010	93,33%	99,99%	99,99%	97,50%				
CPV/FPV Ag CPV + FCoV	2800007	93,33%	99,99%	99,99%	97,50%				
	2800017	93,33%	99,99%	99,99%	97,50%				
	2800011	93,33%	99,99%	99,99%	97,50%	96,00%	97,22%		
CPV + FCoV + Giardia	2800021	93,33%	99,99%	99,99%	97,50%	96,00%	97,22%		
	2800013	95,24%	99,99%						
Dirofilaria / Leishmania	2802015	94,12%	99,99%						
	2802016	94,12%	99,99%	93,62%	97,56%	92,40%	96,67%	95,24%	99,99%
	2802017	94,12%	99,99%	93,62%	97,56%	92,40%	96,67%	95,24%	99,99%
Borrelia	2800022	92,86%	95,83%						
Rotavirus+Bovine Coronavirus+E.coli K99+ Cryptosporidium parvum Ag	2802737	96,15%	98,48%	90,91%	98,77%	90,00%	98,78%	97,92%	97,73%

* Sensitivity and Specificity results from serum testing. Serum is best for testing, whole blood is least accurate. If the results are inconclusive, the test should be repeated with serum.

Physical Examination

High-performance stethoscope with a solid stainless-steel chest-piece for excellent acoustic sensitivity. Features a durable stainless-steel binaural tube with reinforced inner spring, and latex-free Y-type thick-walled tubing with a 3-way sound transmission connector for optimal audio clarity.

Stethoscopes Dual Head Professional

Veterinary Professionals




Item	Description
2801018	● Burgundy
2802392	● Navy Blue

Stethoscopes Dual Head Adult Professional

Item	Description
2802390	● Black
2802391	● Burgundy

Stethoscope Dual Head

Stethoscopes with high density zinc alloy chrome-plated chest-pieces for superior sound transmission. Adjustable chrome-plated brass binaural with white ABS ear tip for optimal comfort.



Item	Description
2801029	● Adult, Black
2801030	● Adult, Grey
2801031	● Adult, Blue
2801032	● Adult, Red
2801033	● Paediatric Small Head
2801034	● Infant Mini Head

Stethoscope Cardiology Professional

Designed for precise auscultation with deep acoustic performance and patient comfort.



Item	Colour
2801024	● Black

Stethoscopes Rappaport Professional

Professional Stethoscope Rappaport/Sprague type with high density zinc alloy chrome-plated chest-pieces for superior sound transmission. Two 560 mm latex-free separately thick wall tubing provides individual sound transmission from chest-piece to both sides of the ear.



Item	Colour
2801020	● Black
2801021	● Navy Blue
2801022	● Neon Orange
2801023	● Royal Blue

Stethoscopes Nurse

Stethoscopes with high density zinc alloy chrome-plated chest-pieces for superior sound transmission. Adjustable chrome-plated brass binaural with white ABS ear tip for optimal comfort.



Item	Colour
2801025	● Black
2801026	● Grey
2801027	● Royal Blue
2801028	● Red

Spare parts



For Dual Head Stethoscope

Item	Description	Quantity
1 2801035	Membrane	1 pc
2 2801036	Eartips Large	1 pair
3 2801037	Eartips Small	1 pair
2 2802393	Ear Tips Large Soft	1 pair
3 2802394	Ear Tips Small Soft	1 pair


For Rappaport Stethoscope

Item	Description	Quantity
1 2801038	Membrane Large	1 pc
2 2801039	Membrane Small	1 pc
3 2801040	Ear Tips Soft	1 pair
4 2801041	Ear Tips Soft Transparent	1 pair

Thermometers

Digital Thermometer

Fast, waterproof thermometer with fever alarm and memory function. Switchable between °C and °F, with replaceable batteries. Ideal for smaller animals.



Item
2800568

Digital Thermometer with Flexi Tip

Designed for veterinary use, this thermometer features a flexible tip. It delivers results in just 10 seconds, is waterproof, and includes a fever alarm. The display is easy to read, with memory function, °C/°F setting, and replaceable batteries.



Item
2800569

Thermometer Cover

Made out of Thin, lubricant-free foil cover that doesn't affect measurement accuracy. Protects against contamination and keeps the thermometer clean. Fits all Covetrus digital thermometers (2800568-2800569).

Item	Quantity
2800832	50 pcs



Percussion Hammers

Percussion Hammer Taylor



Item	Description
2801228	17 cm

Reflex Hammer Troemner



Item	Description
2801242	24 cm

Penlights

Disposable medical diagnostic lights for pupil examination. • Average working time: 3 hours

Item	Description	Quantity
2800793	Penlights Disposable	6 / box



Biopsy Punches

Biopsy Punch, Sterile

The biopsy punches facilitate the removal of skin and tissue samples with the greatest precision and minimal patient trauma. Packed individually.

- The biopsy punches facilitate the removal of skin and tissue samples
- minimal patient trauma
- Packed individually



Item	Size	Quantity
2800971	2.0 mm	10 / box
2800972	3.0 mm	10 / box
2800973	4.0 mm	10 / box
2800974	5.0 mm	10 / box
2800975	6.0 mm	10 / box
2800976	7.0 mm	10 / box
2800977	8.0 mm	10 / box

Dermal Curettes

Dermal Curette, Sterile

The single-sided and extremely sharp cutting edge made from stainless steel enables optimal curettage. The curettes are extremely sharp and are provided with protective caps. It is made from superior stainless steel blade and premium plastic handle. Free from phthalates. Rigid closed-end translucent blade guard protects blade.

- Individually packaged, and ready to use
- Ideal for dermatology, minor surgical procedures, and tissue sampling



Item	Size	Quantity
2801005	2.0 mm	10 / box
2801006	3.0 mm	10 / box
2801007	4.0 mm	10 / box
2801008	5.0 mm	10 / box
2801009	7.0 mm	10 / box

Miscellaneous

Nose Speculum

The flat tapering of the blades provide easier access to tight openings. This nose speculum has an integrated screw joint.



Item	Description
2801227	Hartmann-Halle #2 15.5 cm



Laboratory

Microscope Slides

Microscope Slides Ground Edges

- Washed and polished
- Dimensions: 76 x 26 mm
- Ground edges, 45-degree corner
- Available in plain or twin frosted



Item	Description	Quantity
2800599	Microscope Slides Ground Edges	50 pcs
2802448	With double frosted ends, Ground Edges, Washed and polished	50 pcs

Microscope Cover Glass

Cover slips for microscope slides.

- 18 x 18 mm
- Thickness No.1 (0.13 - 0.16 mm)



Item	Quantity
2800601	100 pcs

Pipette

Pipette Pasteur Polyethylene

Non-sterile with a length of 15 cm and a capacity of 3 ml.

- Capacity: 3 ml
- Length: 15 cm
- Non-sterile, 50 per bag



Item	Quantity
2800002	50 / bag

Applicator

Cotton-Tipped Applicator

Soft, absorbent, cotton-tipped applicators allow for easy application and will not lint. Made from high-quality birch wood. Small: 15 cm, XL: 20 cm.



Item	Description	Quantity
2800165	Wood Small	100 / bag
2800166	Wood XL	50 / bag
2802048	Cotton-Tipped Applicator Double End	200 / bag

Microscope Cleaning Set



Microscope Cleaning Set

Cleaning kit for microscopes, includes: Lens cleaning fluid, Lint free lens tissue, Brush, Air blower, Cotton swabs.

Item	Description
2801953	Microscope Cleaning Set

Patient Care & Animal Handling

Identification

Pet Identification Chip RTU

Sterile pet identification chips with an anti migration coating in a ready to use syringe. Each transponder is individually packed with a disposable applicator.



Item	Size	Country code	Quantity
2800468	Regular (2.12 x 12 mm)	Netherlands	10 / box
2800469	Small (1.5 x 9 mm)	Netherlands	10 / box
2800464	Regular (2.12 x 12 mm)	Germany	10 / box
2800465	Small (1.5 x 9 mm)	Germany	10 / box
2800466	Regular (2.12 x 12 mm)	Austria	10 / box
2800467	Small (1.5 x 9 mm)	Austria	10 / box
2800470	Regular (2.12 x 12 mm)	Poland	10 / box
2800471	Small (1.5 x 9 mm)	Poland	10 / box
2800472	Regular (2.12 x 12 mm)	Universal	10 / box
2800473	Small (1.5 x 9 mm)	Universal	10 / box
2802029	Regular (2.12 x 12 mm)	Denmark	10 / box
2802030	Small (1.5 x 9 mm)	Denmark	10 / box
2802031	Regular (2.12 x 12 mm)	Norway	10 / box
2802032	Small (1.5 x 9 mm)	Norway	10 / box
2802033	Regular (2.12 x 12 mm)	Sweden	10 / box
2802034	Small (1.5 x 9 mm)	Sweden	10 / box

Microchip Reader

Microchip RTU Reader Ring

The chip reader allows you to scan FDX - B (15 digits) and FDX-A (FECAVA) chips used for animal identification according to ISO 11784 standards. Possibility to transfer data and charge the reader via USB cable.

- Internal, ready-to-charge battery
- Scanning time: 8 seconds
- 1000 cycles with one battery
- Dimensions: 135 x 33 mm



Item	Description
2800474	Microchip RTU Reader Ring

Microchip Reader Handheld

- Advanced multi-chip scanner that reads ISO FDX A (not AVID encrypted), FDX B (15 digits), Q5, EM4102, TIRIS HDX and virtually all microchips (transponders) for animal identification according to ISO 11784 standards.
- Dimensions: 155 (L) x 82 (W) x 33 mm (H)
- Weight: 155 g
- Charging mode: mini USB



Item	Description
2800475	MicroChip RTU Reader Handheld
2800476	Battery for Microchip RTU Reader Handheld

Protective Gear

Elizabethan Collar

Elizabethan collar made from soft and flexible polypropylene. Eco-friendly and completely recyclable.



Item	Size	Quantity
2801967	Size 7.5	10 / box
2801968	Size 10	10 / box
2801969	Size 12.5	10 / box
2801970	Size 15	10 / box
2801971	Size 20	10 / box
2801972	Size 25	10 / box
2801973	Size 30	10 / box
2801974	Size 35	10 / box
2801975	Size 40	10 / box

Muzzles

Dog Muzzle for Short-Nosed Breeds

- Adjustable snap-lock clip closure
- Easy to fit securely



Item	Width (cm)	Height (cm)	Size
2802598	20,3	9,8	S
2802599	22,5	11	M
2802600	24,6	13,5	L
2802601	29	15	XL

Protective Boots Petcare Paw

To protect from infection and support the wound healing process.



Size	Breed
S	Yorkshire Terrier, Toy Poodle
M	Jack Russel Terrier, West Highland Terrier
L	Border Collie, Beagle, Dalmatian
XL	Labrador, Boxer, Setter
XXL	Rottweiler, Great Dane, German Shepherd
XXXL	Newfoundlander, St. Bernard, Leonberger

Item	Size
2800821	130 x 60 mm, S
2800822	150 x 75 mm, M
2800823	177 x 85 mm, L
2800824	215 x 90 mm, XL
2800825	230 x 100 mm, XXL
2800826	230 x 110 mm, XXXL

Dog Muzzle Nylon

- Adjustable snap-lock clip closure
- Easy to fit securely



Item	Size
2802305	N.0
2802306	N.1
2802307	N.2
2802308	N.3
2802309	N.4
2802310	N.5
2802311	N.6/4XL
2802312	N.7/5XL
2802313	N.8 Short



Some products from this chapter may only be available in selected countries due to the local registration/regulatory requirements. Please check availability with your Covetrus representative.

Dog Muzzle with leather strap

Item	Size	Circumference	Length
2802615	Size 1	20 cm	6 cm
2802616	Size 2	22 cm	6 cm
2802617	Size 3	24 cm	7 cm
2802618	Size 4	26 cm	7 cm
2802619	Size 5	30 cm	10 cm
2802620	Size 6	31 cm	8 cm
2802621	Size 7	32 cm	10 cm
2802622	Size 8	35 cm	10 cm
2802623	Size 9	38 cm	11 cm
2802624	Size 10	43 cm	14 cm

Dog muzzle with leather strap for Boxers.

Item	Size	Circumference	Length
2800243	Size 4	48 cm	5 cm
2800244	Size 5	56 cm	6 cm



Ectoparasites

Fleacomb



Item	Description
2800665	Without Handle
2800666	With Handle

Tick Remover Pen-Style Plastic



Item	Description
2800667	Pen-Style Tick Remover

Tick Forceps Crossed Steel



Item	Size
2801325	14.5 cm

TICKLESS - Ultrasonic Tick and Flea Repellent


Chemical and fragrance free tick repellent. Works with inaudible soundwaves which irritate ticks and cause them to fall off the pet. Provides at least 6 months of protection and is safe for animals of every age.

Item	Quantity
2802749	1 pc



Nail Clippers

Nail Clippers for different animals



Item	Description	Size
1 2800786	Nail Clipper Guillotine	14 cm
2 2800789	Dog Nail Clipper	12.5 cm
3 2800790	Dog Nail Clipper	16 cm
4 2800787	Cat Nail Clipper	8.5 cm
5 2800788	Cat Nail Clipper	11.5 cm

Miscellaneous

Medical Razor

• Non-sterile



Item	Description	Quantity
2802108	Medical Razor Disposable	100 / box

Bird Ring Scissors



Item	Description	Size
2801054	Bird-ring Scissors	10.5 cm

Underpads

Underpads (Kennel Liners)

Waterproof backing, fast absorbable, non woven, fluff pulp, superabsorbent polymers.



Item	Size	Weight	Quantity
2802502	45 x 60 cm	22 g	30 / pk
2802503	60 x 60 cm	35 g	30 / pk
2802504	60 x 90 cm	55 g	30 / pk

Cadaver Bags

Cadaver Bag Black



Item	Size	Quantity
2802112	XS 67 x 38.5 cm	20 / pk
2802113	S 93 x 57 cm	20 / pk
2802114	M 105 x 77 cm	20 / pk
2802115	L 136.5 x 70.5 cm	20 / pk
2802116	XL 140 x 76 cm	20 / pk

Cat Litter

Cat Litter Bentonite

The cat litter is a natural clay (bentonite). Forms solid and compact lumps when wet.



Item	Description	Size
2800631	Clumping Light-grey	8 kg

Silica Crystals

- Biodegradable
- Antibacterial
- Composition: 100% silicon dioxide
- Ultra-absorbent



Item	Description	Size
2800883	Lumpy Silica	2 x 16 L
2800884	Lumpy Silica	5 L
2800885	Micro Silica	5 L

Recycled Paper Cat Litter

- Made from 100% recycled ink free paper
- Biodegradable
- High odour control
- Exceptionally absorbent



Item	Description	Size
2803129	Recycled Paper Cat Litter	30 L



Large Animal

Cow Ruminal Magnets

Cow Ruminal Magnets

Ruminal magnets prevent hardware disease in cows. The magnets attract ferrous material (nails, wire, screws) in the reticulum and rumen. Designed to sit permanently within the reticulum of the cow: when properly administered, one cow magnet will last the cow's lifetime. The magnets are the ideal size, shape and weight for maximum protection and remedy for bovine traumatic indigestion.



Item	Description	Quantity
1	2802418 Bolus applicator, metal handle	1 pc
2	2802240 Ruminal Cow Magnet Caged	10 / box
3	2800970 Ceramic Ruminal Cow Magnets	12 / box

Marking

Marking Pen/Stick

Crayons (marking pens) based on special waxes and paraffin oil. The marking remains visible on the upper back of cattle and sheep for approx. four weeks, on pigs one to two weeks.



Item	Colour
2800254	Blue
2800255	Green
2800256	Red

Marking Spray

Animal marking spray is a fast-drying and colour-intensive marking dye, suitable for marking cattle and pigs. The marking is visible on pigs for an average of 2 weeks, on cattle somewhat longer.



Item	Colour	Size
2800257	Blue	500 ml
2800258	Green	500 ml
2800259	Red	500 ml

Veterinary Gels

Non-spermicidal Lubricant Gel

- Thick general-purpose lubricating gel
- With skin protection
- Can be used in artificial insemination and obstetrics



Item	Description	Size
1 2800450	Non-spermicidal Gel with special nozzle	500 ml
2 2800451	Non-spermicidal Gel	2.5 L
3 2802926	Non-spermicidal Gel	5 L

Ultrasound Gel

Ultrasound gel for all areas of ultrasound diagnostics and therapy. Bubble-free for exact picturing, skin-friendly, Formaldehyde-free. 5kg is a cubitainer to refill the small dispenser.



Item	Size
1 2802378	250 g
2 2802379	1000 g
3 2802380	5000 g

Teat Instruments

Syringe Mount

Ideal for exact administration of fluid drugs into the teat per syringe.



Item	Description	Quantity
2800260	Luer Non-sterile, 15 mm	100 / box

Milking Tube

Sterile milking tube in blister pack. Luer-lock for easy administration of e.g. antibiotics. Diameter: 3.3 mm.



Item	Description	Quantity
2800785	Luer-lock Sterile, 3,3 mm	20 pcs

Teat Dilator

Teat dilators and teat plugs keep the teat channel open after injuries or operations. They are used for dilatation in case of stenosis.



Item	Description	Quantity
1 2800261	Teat Dilator with Pearl 5.5 cm, Separation Box with Vaseline	20 pcs / box
2 2800262	Teat Dilator in Vaseline 5.0 cm	10 pcs / bag
3 2800263	Teat Plug Wool in Vaseline 5 cm x 6 mm	20 pcs / box



Equipment

Patient Monitoring/Surgery

Thermo Control Professional Heating Mat

Professional heating mat designed to minimise heat loss and support animal recovery during surgery and physiotherapy.

Key features:

- Highly effective FIR infrared heating technology
- Adjustable time and temperature settings (30°C - 42°C)
- Continuous, even heat distribution
- Flexible design suitable for extremities



Item	Description
2802276	With EU cable
2802284	With UK cable
2802302	With AUS cable

Heating Mat Covers

Absorbable microfiber cover with velcro strap.



Item	Description	Quantity
2802285	Heating Mat Cover	3 pcs
2802286	Heating Mat Cover	1 pc

Standard accessories

Item	Description	Quantity
2802556	Controller	1 pc
2802287	Spare Adapter	1 pc

Veterinary Electrocardiograph

- 6.2" LCD touch scree with 1/3 channel rhythm lead display
- Supports real time sampling, pre-sampling, synchronous sampling and sequential sampling
- ECG review, 600 seconds waveforms review and print for any 5 leads within 10 second

Item	Description
2802111	Veterinary Electrocardiograph, EU, UK
2802322	Veterinary Electrocardiograph, AUS

Standard accessories

Item	Description
2802251	5-lead ECG Detachable Clip Connector Veterinary ECG Cable/IEC
2802252	Recording Paper (Roll, 80mm*20m)
2802253	Power Cord (EU) for 2802111
2802553	ECG Softdog Software (incl. Dongle)



Handheld Veterinary Monitor

Handheld Pulse Oximeter & ECG Monitor

The Covetrus Handheld veterinary monitor measures the oxygen saturation of haemoglobin (SpO₂) in arterial blood and the pulse rate. It is equipped with adjustable alarms and also provides visual indicators on a bright 2.4" OLED display. The product also features a gravity activated automatic display rotation for easy viewing in multiple orientations.



Item	Description
2800522	Handheld veterinary monitor with ECG function
2802235	Handheld veterinary pulse oximeter

Standard accessories (included):

Item	Description
1 2800523	Digital SpO ₂ ear clip reusable sensor (9pin VET APK)
2 2800528	3-lead Snap Connector ECG Cable (IEC) (4pin) only for 2800522
3 2802106	Veterinary tongue clip (Size S)
4 2802107	Veterinary tongue clip (Size L)



Optional accessories (not included):

Item	Description
1 2800524	Battery charger stand & AC-DC adaptor for 2800522
2 2800525	Lithium battery for 2800522
3 2800526	Protect silicon jacket (Yellow) for 2800522
4 2800527	Carrying bag (black) for 2800522
5 2800528	3lead Snap Connector ECG cable (IEC, 4PIN)
6 2800529	ECG alligator clip Electrodes (Reusable, Snap Connector, 3pcs/set)



Multi-Parameter Patient Monitor

Veterinary Patient Monitor VPM-25

- 10.4" Color LCD Touch Screen
- 4-12 waveforms
- Measurement Parameter
 - Basic: 3/5 Lead ECG, SpO₂, NIBP, Resp, Temp, HR/PR
 - Option: 12 Lead ECG (Glasgow algorithm), EtcO₂, 2 IBP
- EWS (Early Warning Score): NEWS (National Early Warning Score), NEWS2 and Custom
- OxyCRG
 - Combination of three trend waves, Beat-to-beat heart rate, SpO₂ and Respiration, mainly used in neonatal monitoring
- Communication
 - Infoware G: Intuitive and comprehensive CMS (Central Monitoring System)
 - Infoware X: Highly reliable dedicated communication module to connect via Wi-Fi, WCDMA or LAN



Item	Description
2802733	VPM-25 without Printer
2802734	VPM-25 with Printer

Scope of delivery

Item	Description
2802733/2802734	Veterinary Patient Monitor VPM-25
2801529	Printer Paper (50 mm) *Only When Printer Option is Installed
2801512	SpO ₂ Reusable Sensor WA-103 (Y Sensor)
2801513	SpO ₂ Extension Cable
2801514	Vet Cuff Size 1 (3-6 cm)
2801515	Vet Cuff Size 2 (4-8 cm)
2801516	Vet Cuff Size 3 (6-11 cm)
2801517	Vet Cuff Size 4 (7-13 cm)
2801518	Vet Cuff Size 5 (8-15 cm)
2801519	Vet Cuff Size 6 (17-25 cm)
2802994	ECG 3 Lead Wire (EU Snap Type)
2802974	ECG Cable
2803144	ECG Alligator Clip (3pcs)
2802991	EU Power cord
2801522	UK Power cord
2802992	AU Power cord

Optional accessories (not included):

Capnography Streams



Item	Description
2801509	Side stream LoFlo CO ₂ Module

Power Supply



Item	Description
2802975	1hr battery pack
2802976	4hr battery pack

Temperature Probe



Item	Description
2801524	Temperature probe (rectal)

Straight Sampling Line

These sample lines have a build-in moisture filter.



Item	Description
2801540	Straight sampling line
2801541	Straight sampling line, w/nafion

Oesophageal Probes




Item	Description	Size
2802977	Oesophageal probe	XS
2802978	Oesophageal probe	S
2802979	Oesophageal probe	M

SpO₂ Accessories



Item	Description
2801510	SpO ₂ Rectal Sensor
2801511	SpO ₂ Animal Ear Clip Sensor
2801512	SpO ₂ Animal Sensor (Y-sensor)
2801513	SpO ₂ Extension cable (RCP-058)

Blood Pressure Cuffs



Item	Description	Size	Quantity
2801514	Size 1	3~6cm	20 pcs
2801515	Size 2	4~8cm	20 pcs
2801516	Size 3	6~11cm	20 pcs
2801517	Size 4	7~13cm	20 pcs
2801518	Size 5	8~15cm	20 pcs
2801519	Size 6	17~25cm	1pc


ECG



Item	Description	Quantity
1 2802973	ECG lead wire	
2 2802974	ECG cable	
Item	Description	Quantity
2802315	Veterinary ECG clips	5 pcs
2802109	NIBP Cuff Hose (3m, Rectus) for Vet monitor, XS	1pc

Loflo CO₂ Accessories

Item	Description
2802981	Loflo CO ₂ Nasal Cannula Adult
2802982	Loflo CO ₂ Nasal Cannula Pediatric
2802983	Loflo CO ₂ Nasal Cannula Infant
2802984	Loflo CO ₂ Nasal Cannula w/ O ₂ Adult
2802985	Loflo CO ₂ Nasal Cannula w/ O ₂ Pediatric
2802986	Loflo CO ₂ Nasal Cannula w/ O ₂ Infant
2802987	Loflo CO ₂ Oral/Nasal Cannula Adult
2802988	Loflo CO ₂ Oral/Nasal Cannula Pediatric
2802989	Loflo CO ₂ Oral/Nasal Cannula w/ O ₂ Adult
2902990	Loflo CO ₂ Oral/Nasal Cannula w/ O ₂ Pediatric



Airway Adaptors



Item	Description	Size
2801536	LoFlo Airway adaptor Paediatric	Adult
2801537	LoFlo Airway adaptor Paediatric	Infant
2801538	LoFlo Airway adaptor w/nafion-membrane Paediatric	Adult
2801539	LoFlo Airway adaptor w/nafion-membrane, Paediatric	Infant

Veterinary Glucose Meter

It has been developed specially for dogs and cats and requires samples of at least 0,4 µL from fresh capillary or venous whole blood. The device is designed to simplify the blood glucose measurement.

- Measurement takes only 5 seconds
- Measuring range: 10-600 mg / dl

Kit includes

- Glucose meter
- Carrying case
- 2x batteries
- User manual
- Quick start guide



Item	Description
2800659	VET Glucose Meter mg/dl

Miscellaneous

Tripod on wheels with adjustable height to mount the monitor. The monitor clicks onto the mounting tray and the tripod has a basket to store accessories. With push handle.



Item	Description
2801528	Rolling Stand

Blood Glucose Test Strips

- Avoids an unnecessary long stop for dogs and cats
- Sip-In technique prevents unhygienic blood smearing
- Accurate results
- One type of strips for both glucose meters



Item	Quantity
2800660	25 strips

Handheld Refractometer

Handheld Refractometer for Colostrum

Handheld Brix refractometer which enables quick and reliable evaluation of colostrum quality by estimating IgG levels. With stable IgG concentration between 24-72 hours, the device helps determine whether passive transfer in newborn calves has been successful.

- Measurement range: 0-32% Brix with 0.20% resolution
- Easy evaluation: Brix ≥ 22% indicates high-quality colostrum (>50 mg/ml IgG)



Item	Size
2802110	26 x 40 x 170 mm

Handheld Refractometer for Serum Protein and Urine Specific Gravity

Handheld refractometer for fast and accurate measurement of animal physiological fluids, including urine specific gravity, serum protein and refractive index.

Three-scale measurement:

- serum protein (0-12 g/dl)
- urine specific gravity (1.000-1.040)
- refractive index (1.3330-1.3600)



Item	Size
2802250	26 x 40 x 170 mm

Centrifuges

Digital Centrifuge Benchtop

Table-top digital centrifuge for routine applications such as concentration of cellular suspensions, separation of serum or plasma etc. With a compact and stable design, it presents digital adjustment of centrifugation speed and time.



Item	Description
2802381	Digital Centrifuge Benchtop, EU
2802386	Digital Centrifuge Benchtop, UK
2802388	Digital Centrifuge Benchtop, AUS

Digital Centrifuge Microtubes

Compact and stable digital centrifuge which offers high versatility, easy operation and safety features like overspeed and imbalance detection. With fast acceleration, tool-free rotor changes and precise digital controls via LCD display, it ensures efficient and reliable performance of routine lab applications.



Item	Description
2802382	Digital Centrifuge Microtubes, EU
2802387	Digital Centrifuge Microtubes, UK
2802389	Digital Centrifuge Microtubes, AUS

Accessories



Item	Description
1	2802438 Set of accessories for 2802381, 2802386, 2802388
2	2802439 Carbon brush for centrifuge for 2802381, 2802386, 2802388
3	2802440 Sucking disk type rubber leg for 2802381, 2802386, 2802388
4	2802441 Microhematocrit rotor for 2802382, 2802387, 2802389
5	2802442 Rotor 4 cryotubes 5 ml for 2802382, 2802387, 2802389
6	2802443 PCR rotor, 12 x 1,5/2 ml + 32 x 0,2 ml for 2802382, 2802387, 2802389
7	2802570 Adapter for 13 x 100 mm tubes for 2802381, 2802386, 2802388
8	2802571 Adapter for 13 x 75 mm tubes for 2802381, 2802386, 2802388
9	2802602 PCR 0.2 ml x 12 strips x 4 row for 2802382, 2802387, 2802389
10	2802603 5 ml x 4 tubes for 2802382, 2802387, 2802389

Digital Infusion Devices

Digital Infusion Pump

The multi-function digital infusion pump is easy to use and brings precision, safety and versatility for your day to day fluid therapy requirements. The pump is perfectly adapted for veterinary use with IP34 waterproof design. The 4,3 inch touch screen supports automatic switch between a horizontal and vertical display.

- Uses any pressure infusion set
- Flow rates range from 0,1 ml/h to 1500ml/h with increment of 0,01 ml/h
- Multi-function operation with 4 modes: rate mode, time mode, body weight mode and drip mode
- Foldable handle included

Item	Description
2802674	Digital Infusion Pump, EU
2802675	Digital Infusion Pump, UK
2802676	Digital Infusion Pump, ANZ



Digital Syringe Pump

Our multifunction syringe pump with a 4,3 inch colour touch screen delivers precision, safety, and versatility to support a wide range of veterinary fluid therapy needs. Engineered for clinical reliability, it enables accurate dosing for both small and large animals. Compact and intuitive to operate, it's well-suited for use in both routine treatments and critical care situations.

- Automatic recognition and calibration of syringes (5 ml,10 ml,20 ml,30 ml and 50/60 ml)
- Multi-function operation with 3 modes: rate mode, time mode and body weight mode
- Flow rates range from 0,1 ml/h to 1500 ml/h with increment of 0,01 ml/h
- Stackable design and foldable handle included

Item	Description
2802695	Digital Syringe Pump, EU
2802696	Digital Syringe Pump, UK
2802697	Digital Syringe Pump, ANZ



Digital Infusion Warmer

Heat is transferred via a heat exchanger through the infusion tubing, ensuring continuous warming of the flowing liquid. The target temperature can be precisely set between 30 °C and 50 °C. An intelligent control chip maintains the fluid at the desired temperature, providing stable and consistent warming throughout treatment.

- Compact size and modern appearance
- Easy-to-use for accurate and reliable temperature control

Item	Description
2802671	Digital Infusion Warmer, EU
2802672	Digital Infusion Warmer, UK
2802673	Digital Infusion Warmer, ANZ



Vet Piezo Scaler

Vet Piezo Scaler

- Scaling, endo and automatic water supply functions
- With adjustable intensity and water volume: adjustable in 10 strengths
- Operation with a foot pedal. The device comes with an external 800 ml water tank
- Easy to clean thanks to its smooth surface
- Detachable autoclavable handpiece and tip change:
- Pull out the handpiece and unscrew the scaling tip and endo-chuck after each procedure.
- Pack the handpiece with sterile gauze or sterile bag before sterilizing
- Compatible with external water supply
- Digital controls
- Frequency around 30,000 Hz
- Weight 0.73 kg
- Dimensions 190 x 125 x 60 mm (W x D x H)
- Packaging includes:
- 1 x unit, 1 x clinic kits with 1 x handpiece-E, 1 x tip changer,
- 3 tips (HS-S4-E, HS-P3-E, HS-S1-E)



Item	Description
2801557	Vet Piezo Scaler PS 1
2801831	Vet Piezo Scaler PS 1 LED

Optional Accessories:

Item	Description	Quantity
2801598	Accessory replacement bottle for the Vet Piezo Scaler	1 pc
2801595	S1 scaling tip to remove normal supragingival calculus, also from neck of the tooth	3 pcs
2801596	S4 scaling tip to remove all supragingival calculus, tip design slightly longer than the S1	3 pcs
2801594	P3E perio scaling tip to remove subgingival and deep calculus	3 pcs

Scaler Tips



P3 Perio Tip
P3E perio scaling tip to remove subgingival and deep calculus

Item	Quantity
2801594	3 pcs



S1 Scaling Tip to remove Supragingival Calculus
S1E scaling tip to remove normal supragingival calculus, also from neck of the tooth.

Item	Quantity
2801595	3 pcs



S4 Scaling Tip to remove all Supragingival Calculus
S4E scaling tip to remove all supragingival calculus, tip design slightly longer than the S1.

Item	Quantity
2801596	3 pcs

Handpieces

Optima Turbine with FO Light

The renewed Optima high-speed turbine by BA International features fibreglass light and Smart Coat Technology for improved grip and scratch resistance. With a 20W power output, stainless steel housing, and ceramic bearings, it's suitable for washer disinfectors and autoclaves. It includes a quadruple spray, fibre-optic light, and comes in three fitting options, including KaVo fitting (2801587, BA525LK).



Item	Fitting
2801587	KaVo Fitting 20W Standard Head

Optima Turbine without Light

Fixed base Turbine for Midwest (4 hole) type connections, manufactured by BA International. Made of stainless steel with superior grip.



Item	Fitting
2801589	Midwest 4 Holes 21W Standard Head

Straight Handpieces 1:1

Contra Angle 1:1 without Light/Spray

1:1 contra angle handpiece with push-button bur release and standard E-fitting. It features a smooth ergonomic design for better grip and easy cleaning. The handpiece is autoclavable, thermally disinfectable, and operates without light or spray.



Item	Description
2801571	BA45E - Push Button

Optima Straight 1:1 Without Spray

Straight handpiece by BA with a 1:1 direct ratio, no spray or light. Made with a standard E-fitting and suitable for washer disinfectors and autoclaves.



Item	Description
2801576	BA175E without Spray

Contra Angle 1:1 Latch Type without Light

Ultimate Contra Angle with 1:1 direct ratio and latch-type bur release. Operates without spray or light, features a standard E-fitting, and is autoclavable.



Item	Description
2801568	BA101 - Latch type

Ultimate Contra Angle Prophy 1:1

Ultimate Contra Angle handpiece for prophy cups and brushes, available in Snap-on type (2801570). It features a 1:1 direct ratio, standard E-fitting, and is autoclavable.



Item	Description
2801570	BA 103 for Snap-on Cups & Brushes

Water Attachments

Water Attachment for Contra Angle Handpieces

Water clip for BA100 series contra-angle handpieces.

- Fits:
- BA101 (2801568)
 - BA103 (2801570)

Item
2801577



LED Curing Light

Cordless, lightweight pen-style curing light with an output of over 1200 mW/cm², ideal for curing light-curing composites. It features a digital control panel with an audible signal for the last 3 seconds and is easy to store and handle. The high-performance diode cures up to 3 mm of resin-based material in 10 seconds. Offers three curing modes - continuous, ramp-up, and pulse - with adjustable timing. Operates at 230 V, 50/60 Hz, with a wavelength of 450-470 nm and a long-life 5W blue LED. The handpiece weighs 154 g and measures 25 x 3.5 x 2.5 cm; charger unit is 11.5 x 8 x 7.5 cm. Includes an 8 mm fiber-optic guide, two small orange tips, an orange shield, and a table charger.



Item **Description**
2801597 Curing Light LED 1200 Cordless Silver

Thermosealer

Electronic thermosealing machine for continuous cycle processing. The machine is using the latest generation of heating elements with the possibility to adjust the sealing temperature. For use with sterilisation rolls.

- Compact design with sealing lever positioned at the front
- The machine creates a 12 mm sealing strip and has a sealing are width of 310 mm
- Has an integrated cutter



Item **Description**
2801565 Thermosealer BA Ultiseal BASE740

Ultrasonic Cleaner

Compact ultrasonic scaler unit specifically designed for veterinary dentistry.

- Scaling, endo and automatic water supply functions
- With adjustable intensity and water volume: adjustable in 10 strengths
- Operation with a foot pedal. The device comes with an external 800 ml water tank
- Detachable autoclavable handpiece and tip change: Pull out the handpiece and unscrew the scaling tip and endo-chuck after each procedure. Pack the handpiece with sterile gauze or sterile bag before sterilizing



Item **Description**
2801559 BA Ulticlean Ultrasonic Bath (UK)
2801560 BA Ulticlean Ultrasonic Bath (EU)
2801564 Full Basket for Ultrasonic bath (size 25 x 13 x 7.3 cm)

Microscopes

Microscope Binocular

Microscope for routine analysis. Can be used for simple cytology and histopathology, urine, blood, smears, parasites, micro-organisms, skin scrapings, ear mite and fungus.



Item **Description**
2801854 Microscope Binocular
2801887 Aluminium case for Microscope Binocular

Microscope Trinocular

Microscope for research and laboratories. Can be used for cytology and basic histopathology, urine, blood, smears, parasites, micro-organisms, skin scrapings, ear mite and fungus.



Item **Description**
2801855 Microscope Trinocular
2801888 Aluminium case for Microscope Trinocular

Microscope Trinocular HE

Microscope with state-of-the-art features, a high resolving power of the objectives and the 3W NeoLed transmitted LED illumination designed for all kind of demanding applications as well for fundamental research.



Item **Description**
2802128 Microscope Trinocular HE
2801888 Aluminium case for Microscope Trinocular

Microscope WiFi Camera

Adding a WiFi camera to your microscope is ideal for wireless sharing/viewing microscope images.



Item **Description**
2802932 Microscope Wifi Camera

Laboratory Microscope Cleaning Set

Item
2801953

PAGE
59

Microscope Comparison Overview

Feature	Binocular 2801854	Trinocular 2801855	Trinocular HE 2802128
Optical System	Achromatic semi-plan	Infinity corrected optical system (IOS)	Infinity corrected optical system (IOS)
Head Type	Binocular head, 30° inclined, 360° rotatable	Siedentopf trinocular, 30° inclined, IPD 48-76 mm, diopter adjustable, 2 height positions, 360° rotatable	Siedentopf trinocular, 30° inclined, IPD 47-76 mm, diopter adjustable
Eyepieces	WF 10x/20 (with eye guard, anti-mold)	WF 10x/20 (with eye guard, anti-mold)	EWF 10x/22 (anti-mold)
Nosepiece	Quadruple (inward)	Quintuple	Quintuple
Objectives	Semi-Plan Achromatic: 4x/0.1, 10x/0.25, 40x/0.65 (spring), 100x/1.25 (oil, spring)	E-Plan Achromatic IOS: 4x/0.1, 10x/0.25, 40x/0.65 (spring), 100x/1.25 (oil, spring)	E-Plan IOS: 4x/0.1, 10x/0.25, 40x/0.65 (spring), 100x/1.25 (oil, spring)
Stage	Mechanical stage 130 x 130 mm, travel range 70 x 28 mm	Rackless stage 197/152 x 131 mm, travel range 75 x 36 mm	Rackless stage 187/230 x 140 mm, double slide holder, integrated X-Y stage 79 x 52 mm
Focusing	Coaxial coarse/fine, 20 mm range, 2 µm graduation, 0.4 mm/rotation	Coaxial coarse/fine, 20 mm range, 1.11 µm precision, 200 µm/rotation	Coaxial coarse/fine, 20 mm range, 1 µm graduation, 200 µm/rotation
Illumination	1 W NEOLED™	3 W NEOLED™	3 W NEOLED™
Dimension (HxWxD)	360 x 176 x 227 mm	409 x 195 x 235 mm	407 x 210 x 255 mm
Weight	4.0 kg	6.5 kg	9.0 kg
Warranty	10 years warranty on mechanical parts and optics (2 years on electronical parts)		

Scales

Compact Scale

This Compact Scale is ideal for small pets such as kittens, rabbits, rodents, birds, and reptiles.

- Dimensions: 19.5 cm x 13 cm x 2.8 cm
- Max. 3 kg, Division 1g
- Clear LCD screen, Tare function



Item	Description
2802634	Compact Scale for Lab or small animals

Pet Scale

This Pet Scale is perfect for weighing smaller dogs, cats, and other similar-sized animals. Its wide, low-profile surface ensures stability and comfort for pets during weighing. Despite its size, the scale is remarkably lightweight, making it effortless to store or move between locations.

- Dimensions 70 cm x 32 cm x 4.25 cm
- Max. 30 kg, Divisions 2 g
- Clear LCD screen, Tare function



Item	Description
2802635	Pet Scale for smaller dogs and cats

Clippers and Trimmers

Premium Pet Hair Trimmer

The trimmer for the specialist areas such as paws face and for precision contouring.

- Very good cutting results due to clipping head made from carbon steel
- Works very quietly and without vibration
- Low weight of just 115 g incl. battery and blade



Item	Description
1 2802506	Premium Trimmer
2 2802507	Replacement blade

Pet Hair Clipper

Compact and powerful pet hair clipper with safe & sharp clipping blade.




Item	Description
1 2802383	Pet Hair Clipper
2 2802384	Comb small 3/6mm
3 2802385	Comb medium 9/12mm
4 2802395	Standard size blade
5 2802396	Narrow size blade
6 2802397	USB charging cable



Instruments

Needle Holders

Needle Holder Mathieu



Item	Description	Size
1 2801201	Needle Holder Mathieu	14 cm
2 2801200	Needle Holder Mathieu	17 cm
3 2801202	Needle Holder Mathieu Slim	14 cm
4 2801203	Needle Holder Mathieu Strong	17 cm
5 2801204	Needle Holder Mathieu Strong	20 cm

Needle Holder Olsen-Hegar



Item	Size
2801218	21 cm
2801217	19 cm
2801221	18 cm
2801216	17 cm
2801222	14 cm

Needle Holder Baumgartner



Item	Size
2801190	14 cm

Needle Holder Mayo-Hegar



Item	Size
2801208	14 cm
2801205	16 cm
2801209	18 cm
2801210	20 cm
2801941	24 cm

Needle Holder Micro Castroviejo



Item	Description	Size
2801211	with Lock Straight	14 cm
2801212	with Lock Straight	16 cm



Item	Description	Size
2801213	with Lock Straight Ø 8 mm	12 cm

Needle Holder Brown



Item	Size
2801191	13 cm

Needle Holder Crile-Wood



Item	Size
2801193	15 cm
2801194	18 cm
2801195	20 cm

Needle Holder Derf



Item	Size
2801196	12 cm

Needle Holder Micro Castroviejo TC



Item	Description	Size
2801214	With Lock Straight Serrated TC	14 cm

Needle Holder Crile-Wood TC



Item	Size
2801192	15 cm

Needle Holder Gillies Left Handed



Item	Size
2801197	16 cm

Needle Holder Gillies Right Handed



Item	Size
2801198	16 cm

Premium Range with TC Insert



Needle Holder Mayo-Hegar TC



Item	Size
2801206	16 cm
2801207	20 cm
2801940	24 cm

Needle Holder Ryder-Vascular TC



Item	Size
2801223	13 cm
2801224	15 cm
2801225	18 cm

Needle Holder Olsen-Hegar TC



Item	Size
2801219	14 cm
2801220	16.5 cm

Scalpel Handles & Blades

Scalpel Handle #3



Item	Description	Size
2801253	Suitable for 1 Blade	3

Scalpel Handle #4



Item	Description	Size
2801254	Fits all larger sizes of scalpel blades	4

Scalpel Handle #3 Round Handle

Size 3 fits scalpel blades #10-19.

Item	Description	Size
2801050	Scalpel Handle #3 Round Handle	3



Disposable Scalpels

Sterile disposable scalpels made of stainless steel with plastic handle. They are characterised by extraordinary stability, exact cutting edges and long-lasting degree of shaft.

Item	Size	Quantity
1 2800978	#10	10 pcs
2 2800979	#11	10 pcs
3 2800980	#12	10 pcs
4 2800981	#15	10 pcs
5 2800982	#20	10 pcs
6 2800983	#21	10 pcs
7 2800984	#22	10 pcs
8 2800985	#23	10 pcs
9 2800986	#24	10 pcs



Scalpel Blades



Scalpel Blades

High-quality scalpel blades are manufactured using tightly controlled sharpening and heat treatment processes to ensure exceptional sharpness, strength, and durability.

Each blade features a precisely honed cutting edge that maintains its sharpness over extended use, ensuring clean and accurate incisions. The blades conform to British Standard BS 2982:1992 and ISO 7153-1:1991, guaranteeing consistent quality and safety.

- Individually packaged in sterile aluminium foil for maximum hygiene and convenience
- Additionally protected with VPI (Volatile Corrosion Inhibitor) paper
- Blade size clearly marked on the reverse side of each peel-pack

	Carbon Steel	Stainless Steel	Description	Quantity
1	2800987	2800996	#10	100 pcs
2	2800988	2800997	#11	100 pcs
3	2800989	2800998	#12	100 pcs
4	2800990	2800999	#15	100 pcs
5	2800991	2801000	#20	100 pcs
6	2800992	2801001	#21	100 pcs
7	2800993	2801002	#22	100 pcs
8	2800994	2801003	#23	100 pcs
9	2800995	2801004	#24	100 pcs

Stitch Cutters Carbon Steel

Precision sharpening and heat treatment ensure consistently sharp, durable blades. Rounded blade tips enhance patient comfort and reduce percutaneous injury risk. Each blade is individually packed in aluminum foil with VPI paper, supplied sterile in unit boxes of 100. Complies with British Standard 2982:1992, has a 5-year shelf life, and is gamma sterilized.



Item	Length	Quantity
2801542	9.3 cm	100 pcs

Long Stitch Cutters Stainless Steel

The long, pencil-like stitch cutter features a classic ribbed grip for better security than the thumb version. Its ergonomic handle offers stability and control for both left- and right-handed users, while ensuring a clear view of the surgical site. Precision sharpening and heat treatment deliver consistently sharp, durable blades. It is phthalate-free and not made with natural rubber latex. A rigid, closed-end blade guard helps prevent sharps injuries. Do not recap the guard.



Item	Description	Quantity
2801543	Steel with Handle	10 pcs

Scissors

Explanation of # numbers in product descriptions



Scissors Curved #1



Item	Size
2801261	14 cm
2801262	16 cm
2801263	18 cm

Scissors Straight #1



Item	Size
2801294	11 cm
2801295	14 cm
2801296	16 cm
2801297	18 cm

Scissors Curved #2



Item	Size
2801264	14 cm

Scissors Straight #2



Item	Size
2801298	13 cm
2801299	14 cm
2801300	16 cm

Scissors Curved #3



Item	Size
2801265	14 cm

Scissors Straight #3



Item	Size
2801301	13 cm
2801302	14 cm
2801303	16 cm

Surgical Scissors Curved



Item	Description	Size
2801313	Surgical Scissors Curved #1	13 cm



Item	Description	Size
2801316	Surgical Scissors Curved #2	13 cm
2801304	Surgical Scissors Curved #2	14 cm



Item	Description	Size
2801314	Surgical Scissors Curved #3	13 cm

Surgical Scissors Straight



Item	Description	Size
2801318	Surgical Scissors Straight #1	13 cm



Item	Description	Size
2801317	Surgical Scissors Straight #2	11 cm
2801305	Surgical Scissors Straight #2	14 cm



Item	Description	Size
2801319	Surgical Scissors Straight #3	11 cm

Ophthalmic Scissors

Iris Curved



Item	Description	Size
2801271	Scissors Iris Curved #3	9 cm
2801315	Scissors Iris Curved #3	10 cm
2801267	Scissors Iris Standard Curved #3	11 cm

Iris Straight



Item	Description	Size
2801273	Scissors Iris Standard Straight #3	11 cm
2801274	Scissors Iris Standard Straight #3	10 cm
2801275	Scissors Iris Straight #3	9 cm

Iris Standard Curved

Item	Description	Size
2801268	Scissors Iris Standard Curved #3	9 cm
2801270	Scissors Iris Standard Curved #3	10 cm



Specialised Scissors

Mayo Scissors



Item	Description	Size
2801281	Scissors Mayo Curved	14 cm
2801280	Scissors Mayo Curved	17 cm
2801283	Scissors Mayo Straight	14 cm
2801282	Scissors Mayo Straight	17 cm

Metzenbaum Scissors



Item	Description	Size
2801291	Scissors Metzenbaum Curved #1	14 cm
2801288	Scissors Metzenbaum Curved #1	18 cm
2801290	Scissors Metzenbaum Straight #1	14 cm
2801289	Scissors Metzenbaum Straight #1	18 cm

Mayo-Stille

Item	Description	Size
2801284	Scissors Mayo-Stille Curved	15 cm
2801285	Scissors Mayo-Stille Curved	17 cm
2801286	Scissors Mayo-Stille Curved	19 cm



Other Scissors

Goldmann-Fox Curved Toothed



Item	Description	Size
2801266	Scissors Goldmann-Fox Curved Toothed	13 cm

Stevens Straight



Item	Description	Size
2801292	Scissors Stevens Straight #1	10 cm

Kelly Straight



Item	Description	Size
22801277	Scissors Kelly Straight	16 cm

Ligature Spencer Straight



Item	Description	Size
2801278	Scissors Ligature Spencer Straight	12.5 cm

Bird-ring Scissors

Suitable for aluminium ring removal.



Item	Description	Size
2801054	Bird-ring Scissors	10.5 cm

Stevens Straight



Item	Description	Size
2801293	Scissors Stevens Straight #3	10 cm

Suture Scissors Littauer Straight



Item	Description	Size
2801320	Suture Scissors Littauer Straight	14 cm

Suture Scissors Spencer Straight



Item	Description	Size
2801321	Suture Scissors Spencer Straight	9 cm
2801322	Suture Scissors Spencer Straight	10.5 cm

Nail Scissors Curved



Item	Description	Size
2801188	Nail Scissors Curved	10.5 cm

Scissors All Purpose Plastic Handle Black



Item	Size
2801255	15 cm
2801256	18 cm

Bandage Scissors



Item	Description	Size
2801257	Scissors Bandage Lister	9 cm
2801258	Scissors Bandage Lister	11 cm
2801259	Scissors Bandage Lister	14 cm
2801260	Scissors Bandage Lister	18 cm
2801279	Scissors Bandage Lister with Pocket Clip	14 cm

Supercut Scissors and TC Scissors

Iris Supercut



Item	Description	Size
2801272	Scissors Iris Curved Supercut	11 cm
2801276	Scissors Iris Straight Supercut	11 cm

Mayo-stille Curved TC



Item	Description	Size
2801287	Scissors Mayo-stille Curved TC	17 cm

Forceps

Hemostat Forceps

Rochester Pean

Curved



Item	Description	Size
2801152	Curved	14 cm
2801153	Curved	16 cm
2801154	Curved	18 cm
2801149	Curved	20 cm

Straight



Item	Description	Size
2801156	Straight	14 cm
2801157	Straight	16 cm
2801150	Straight	18 cm

**Rochester-Pean
Delicate Curved**



Item	Description	Size
2801151	Hemostat Forceps, Delicate Curved	13 cm

Delicate Straight



Item	Description	Size
2801155	Hemostat Forceps, Delicate Straight	13 cm

**Kelly
Curved**



Item	Description	Size
2801142	Hemostat Forceps Kelly Curved	14 cm

Straight



Item	Description	Size
2801143	Hemostat Forceps Kelly Straight	14 cm

**Kocher
Delicate Curved**



Item	Description	Size
2801145	Hemostat Forceps Kocher Delicate	13 cm
2801144	Hemostat Forceps Kocher Delicate	14 cm

Delicate Straight



Item	Description	Size
2801146	Hemostat Forceps Kocher Delicate	13 cm
2801147	Hemostat Forceps Kocher Delicate	14 cm
2801148	Hemostat Forceps Kocher Delicate	16 cm

**Ochsner-Kocher
Straight**



Item	Description	Size
2801125	Forceps Ochsner-Kocher Straight	14 cm
2801126	Forceps Ochsner-Kocher Straight	16 cm
2801123	Forceps Ochsner-Kocher Straight	18 cm

Curved



Item	Description	Size
2801124	Forceps Ochsner-Kocher Curved	18 cm

**Halstead Mosquito
Curved**



Item	Description	Size
2801141	Hemostat Forceps Halstead	12.5 cm

Straight



Item	Description	Size
2801139	Hemostat Forceps Halstead	12.5 cm

Straight 1 x 2 Teeth



Item	Description	Size
2801140	Hemostat Forceps Halstead Mosquito Straight 1 x 2 Teeth	12 cm

**Mosquito Forceps
Curved**



Item	Description	Size
2801179	Mosquito Forceps Curved	10 cm

Straight



Item	Description	Size
2801180	Mosquito Forceps Straight	10 cm

**Crile
Curved**



Item	Description	Size
2801137	Hemostat Forceps Crile Curved	14 cm

Straight



Item	Description	Size
2801138	Hemostat Forceps Crile Straight	14 cm

Hartmann Forceps

Curved



Item	Description	Size
2801121	Forceps Hartmann Curved	10 cm

Straight



Item	Description	Size
2801122	Forceps Hartmann Straight	10 cm

Spencer-Wells Forceps

Curved



Item	Description	Size
2801127	Forceps Spencer-Wells Curved	14 cm
2801129	Forceps Spencer-Wells Curved	18 cm

Straight



Item	Description	Size
2801128	Forceps Spencer-Wells Straight	14 cm

Tissue Forceps

Tissue Forceps Allis



Item	Description	Size
2801338	4 x 5 Straight	12.5 cm
2801339	5 x 6 Straight	15 cm
2801340	5 x 6 Straight	19 cm

Tissue Forceps 1 x 2 Teeth



Item	Size
2801327	13 cm
2801328	14.5 cm
2801343	16 cm
2801326	18 cm

Tissue Forceps Straight 1 x



Item	Description	Size
2801342	Tissue Forceps Straight 1 x 2	14 cm

2 x 3 Teeth



Item	Description	Size
2801329	Tissue Forceps 2 x 3 Teeth	11 cm
2801344	Tissue Forceps 2 x 3 Teeth	14 cm

Tissue Forceps Adson-Brown Surgical 7 x 7



Item	Description
2801336	12 cm Surgical 7 x 7
2801330	15 cm 7 x 7 Teeth

Other Tissue Forceps

Adson-Brown 9 x 9 Teeth



Item	Description	Size
2801331	Tissue Forceps Adson-Brown, 9 x 9 Teeth	12 cm

Treves Forceps Surgical 1x2



Item	Description	Size
2801132	Forceps Treves Surgical 1 x 2	12.5 cm

Tissue Forceps Straight



Item	Description	Size
2801345	Tissue Forceps Straight	18 cm

Anatomical



Item	Description	Size
2801120	Forceps Anatomical Serrated	11 cm
2801131	Forceps Anatomical Straight	14 cm

Adson-Brown



Item	Description	Size
2801332	Tissue Forceps Adson 1 x 2 Surgical	12 cm

Adson Mini



Item	Description	Size
2801333	Tissue Forceps Adson Mini	12 cm

Adson Surgical 1 x 2



Item	Description	Size
2801334	Tissue Forceps Adson Surgical 1 x 2, Smooth with Tying Platform	13 cm

Adson Surgical Fine Micro



Item	Description	Size
2801335	Tissue Forceps Adson Surgical, Fine Micro	12 cm

Dressing Forceps

Straight



Item	Size
2801079	13 cm
2801078	14.5 cm
2801080	16 cm
2801081	20 cm

Adson-Micro Anatomical



Item	Description	Size
2801337	Tissue Forceps Adson-Micro, Anatomical	12 cm

Graefe 1 x 2




Item	Description	Size
2801341	Tissue Forceps Graefe 1 x 2	10 cm

Adson Anatomical



Item	Description	Size
2801077	Dressing Forceps Adson Anatomical	12 cm

Eye Forceps/Ophthalmic Forceps



Eye Forceps

Item	Description	Size
1 2801115	Dressing Iris, Straight	10.5 cm
2 2801114	Dressing Iris Straight	10.5 cm
3 2801119	Serrated, Curved	10.5 cm



Graefe Straight 1.0 mm

Item	Description	Size
2801118	Eye Forceps Graefe Straight 1.0 mm	10 cm



Graefe

Item	Description	Size
1 2801116	Eye Forceps Graefe with Latching	10 cm
2 2801117	Eye Forceps Graefe without Latching	11 cm

Intestinal Forceps

Doyen Straight

Item	Description	Size
2801171	Intestinal Forceps Doyen Straight	23.5 cm



Bone Rongeur



Item	Description	Size
2801246	Rongeur forceps Luer slim curved	15 cm

Ear Forceps



Item	Description	Size
2801082	Ear Forceps Troeltsch Anatomical	12 cm

Swab Forceps



Item	Description	Size
2801323	Foerster-Ballenger Straight	20 cm

Tongue Forceps




Item	Description	Size
2801346	Tongue Forceps Young	16 cm

Splinter Forceps Feilchenfeld



Item	Description	Size
2801308	Splinter Forceps Feilchenfeld	9 cm
2801309	Splinter Forceps Feilchenfeld	13.5 cm

Splinter Forceps Stieglitz



Item	Description	Size
2801310	Curved	13.5 cm
2801311	Straight	13.5 cm

Pliers

Plier College



Item	Description	Size
2801238	Plier College #3	15 cm

Plier Meriam



Item	Description	Size
2801239	Plier Meriam #1	16 cm

Ear Plier



Item	Description	Size
2801083	Ear Plier Lucae Anatomical	14 cm

Grasping Forceps

Hartmann Serrated

Item	Size
1 2801136	8.5 cm
2801133	14 cm
2801134	16 cm
2 2801135	1 x 2



Splinter Forceps Straight



Item	Size
2801130	9 cm

Splinter Forceps Straight Serrated



Item	Size
2801312	11 cm

Towel Clamps

Towel Clamp Backhaus



Item	Size
2801353	8 cm
2801354	11 cm
2801355	13 cm
2801352	15 cm

Towel Clamp Jones



Item	Size
2801356	5.5 cm

Towel Clamp Schaedel Straight



Item	Size
2801357	9 cm

Wound Retractors

Wound Retractor Gelpi-Loktite



Item	Size
2802019	9 cm
2801360	13.5 cm

Wound Retractor Weitlaner Blunt



Item	Size
2801361	13.5 cm

Retractor Farabeuf



Item	Description	Size
2801243	Retractor Farabeuf #1 #2	11.5 cm

Retractor Senn



Item	Description	Size
2801244	Blunt 22 x 7 mm	16 cm
2801245	Sharp 22 x 7 mm	16 cm

Bone Instruments

Bone Curette Martini



Item	Description	Size
2801058	Bone Curette Martini #0/1	14 cm

Bone Curette Willinger



Item	Description	Size
2801061	Bone Curette Willinger Oval #0/0	17 cm
2801062	Bone curette Willinger Oval #1/1	17 cm

Bone Curette Volkmann



Item	Description	Size
2801059	Bone Curette Volkmann #0	17 cm



Item	Description	Size
2801060	Bone Curette Volkmann #0000	17 cm

Bone Curette Hemingway



Item	Size
2801063	17 cm

Bone Rongeur Forceps Mini Friedman



Item	Size
2801065	12.5 cm

Other Products

Percussion Hammer Taylor



Item	Size
2801228	6 3/4"

Reflex Hammer Troemner



Item	Size
2801242	24 cm

Spay Hook



Item	Description	Size
2801066	for Cats	20 cm
2801899	for Small-Medium Dogs	19 cm

Nose Specula



Item	Description	Size
2801227	Hartmann-Halle #2	15.5 cm

Needle Deschamps



Item	Description	Size
2801189	Right Handed	20 cm

Instruments for Large Animals

Birthing-chain

Used for aiding in the delivery of a foetus.



Item	Size
2801056	150 cm

Birthing-handle

Complementary with Birthing-chain (2801056).



Item	Size
2801057	13 cm

Hoof Rasp

Item	Description	Size	Weight
2801159	with wooden handle	52 cm	750 g



Rough surface



Fine surface

Hoof Examination Forceps



Item	Size	Weight
2801158	48 cm	1190 g

Claw Scissors for Sheep



Item	Size
2801068	23 cm

Instrument Boxes & Accessories

Instrument Box - Perforated

Made from aluminium, autoclavable and lightweight. Both lid and bottom are perforated.



Item	Description	Size
2801937	Small	300 x 140 x 40 mm
2801938	Medium	300 x 140 x 70 mm
2801939	Large	300 x 140 x 100 mm

Permanent Filter

Permanent filters are made from Teflon. They are suitable for our perforated instrument boxes (2801937, 2801938 and 2801939) and the spare lid (2801943). Suitable for 1200 sterilisation cycles.



Item	Description	Size
2802071	For Perforated Instrument Boxes	235 x 118 mm

Indicator Labels

Indicator labels suitable for perforated instrument boxes (2801937, 2801938 and 2801939) and the spare lid (2801943). Made from paper, disposable. Packed in boxes of 100 pieces.



Item	Description	Size
2802072	For Perforated Instrument Boxes	100 pcs

Kidney Dish



Item	Size
2801173	170 x 100 x 35 mm
2801174	250 x 140 x 40 mm
2801175	275 x 150 x 45 mm
2801176	200 x 100 x 40 mm

Instrument Box with Lid - Non-Perforated



Item	Description	Size
2801160	Without Holder on Lid	18 x 12 x 5 cm
2801161	With Holder on Lid	20 x 10 x 5 cm
2801163	With Holder on Lid	22 x 15 x 5 cm
2801164	Without Holder on Lid	26 x 15 x 5 cm
2801165	Without Holder on Lid	30 x 20 x 5 cm

Silicone Mat

Secure storage for instruments.



Item	Size
2801942	270 x 125 mm

Instrument Tray



Item	Size
2801162	21 x 16 x 1 cm

Instrument Cassette



Item	Description	Size
2801169	5 Instruments including 3 Silicone Rails	18 x 7 x 2.5 cm
2801167	10 Instruments including 3 Silicone Rails	18 x 14.5 x 2.5 cm
2801168	16 Instruments including 6 Silicone Rails	28 x 18 x 3.5 cm

Instrument Tapes

The marking tapes are for easy labelling and identification of instruments. Autoclavable and suitable for washer disinfection.

The kit includes:

- 2 x rolls Neon Purple
- 2 x rolls Neon Pink
- 1 x roll Neon Blue
- 1 x roll Neon Yellow
- 1 x roll Neon Green
- 1 x roll Neon Orange
- Organiser Rack



Item	Description	Size
2801858	ID Tape Kit Assorted Neon Colours	8 x 1 m
2801859	● ID Tape Neon Blue	1 x 3 m
2801860	● ID Tape Neon Green	1 x 3 m
2801861	● ID Tape Red	1 x 3 m
2801862	○ ID Tape White	1 x 3 m
2801863	● ID Tape Neon Yellow	1 x 3 m
2801864	● ID Tape Grey	1 x 3 m
2801865	● ID Tape Copper	1 x 3 m
2801914	● ID Tape Neon Orange	1 x 3 m
2801915	● ID Tape Neon Purple	1 x 3 m
2801916	● ID Tape Neon Pink	1 x 3 m

Instrument Rings

Instrument ID rings in size: 1/4" (0.64 cm) and neon (vibrant) colours. Autoclavable up to 36°C/275°F. Suitable for washer-disinfection. Not made with natural rubber latex.

The set includes:

- 2 x 25 Neon Purple
- 2 x 25 Neon Pink
- 1 x 25 Neon Blue
- 1 x 25 Neon Yellow
- 1 x 25 Neon Green
- 1 x 25 Neon Orange
- Dispensers and rack to keep them organised and accessible for use



Item	Description	Size
2801857	ID Rings Assorted, Neon Colours	8 x 25 pk Rings

Maintenance

Oil Spray

Universal oil spray for straight and contra-angle handpieces. Degreases, cleans, and lubricates. Including spray head.



Item	Size
2801492	500 ml



Dentistry

Hygiene & Maintenance

Dental Unit & Scaler Solution

Dental solution with chlorhexidine that serves as a coolant for the moving parts of turbines, for handpieces and contra-angles. It reduces the cutting time of veterinary dental instruments. Can also be used as an odourant against bad breath (flush in the mouth to get fresh breath).



Item	Size
2801497	250 ml

Sharpening Stone

Wooden box not included.



Item	Description	Size
2801306	Arkansas Fine Grit Rectangular	15 x 5 cm

Dental Care

Dental Care Kit

The veterinary Dental Care kit consists of a red finger toothbrush for regular brushing and a blue finger toothbrush for massaging the gum.



Item	Quantity
2800664	2 pcs

Dental Cleaning Kit

The Veterinary Dental Cleaning kit contains 60g toothpaste and a toothbrush.



Item	Description
2800818	Toothpaste + Toothbrush

Fluoride Foam Vanilla

Topical oral fluoride foam with 1.5% fluoride and vanilla flavour, gluten free. Apply on clean teeth.



Item	Net weight
2800812	125 g

Fluoride Gel

Berry-flavored gel with 0.4% stannous fluoride for at-home pet dental care. Apply a pea-sized amount twice weekly with a soft or finger toothbrush.



Item	Net weight
2800813	120 g

Dental Chlorhexidine Paste

0.12% Chlorhexidine gel with thick formula for easy, lasting application. Use pea-sized amount once or twice daily. Ideal for short-term home care after cleaning; suitable as retail product.



Item	Net weight
2800815	114 g

Dental Chlorhexidine Spray

Oral Chlorhexidine Spray with 0.12% Chlorhexidine and alcohol free. It can be used once or twice a day for oral hygiene maintenance.



Item	Size
2800816	30 ml

Dental Instruments

Mouth Retractors

Mouth retractors for dogs and cats with spring and self-locking.

- Easy to use
- Plastic/rubber inserts



Item	Description	Size
2801183	For Cats	6 cm
2801186	For Dogs, Small	7 cm
2801185	For Dogs, Medium	10.5 cm
2801184	For Dogs, Large	15 cm

Mouth Mirror Handle - Stainless Steel



Item	Description	Size
2801181	Mouth Mirror without Handle	14 cm
2801182	Mouth Mirror Handle	14 cm

Mouth Mirror Fahrenkrug



Item	Description
2801351	For Horses and Large Animals

Dental Explorer

Double-ended explorer used to test tooth stability and the condition of the gum margins.



Item	Description	Size
2801076	#17/23 Double-ended Comfort	16.5 cm

Tooth Probe

Tooth probes have a thin blunted end and are used to look for caries and calculus detection. With octagonal handle and thin unbreakable tip. The tooth probes have a single-ended handle solid.

Item	Description
1 2801348	#9 Single-ended Octagonal Handle Solid Serrated
2 2801349	#TU17 Wilkins/Tufts Single-ended Octagonal Handle Solid Serrated



Periodontal Probe

A periodontal probe is a dental instrument used in the diagnosis of gum diseases. They have a thin blunted end and markings at the tip. They are used to measure the periodontal pocket depth (space between the gums and the tooth).

Item	Description
1 2801229	#CP11 Graduate 3-6-8-11 mm Comfort
2 2801231	#CP12 Graduate 3-6-9-12 mm Comfort
3 2801232	#CP15 North Carolina Graduate 1-2...-14-15 mm Comfort
4 2801233	#CP2 Graduate 2-4-6-8-10-12 mm Comfort



Hand Scaler

Scalers are useful for the removal of supragingival calculus and plaque (above the gum line).

Item	Description
1 2801248	#H6-H7 Hygienist Double-ended Comfort
2 2801249	#CK6 Crane-Kaplan Double-ended Comfort
3 2801250	Jacquette
4 2801251	With Wings Left 15 cm
5 2801252	With Wings Right 15 cm



Gracey Curette, Double-ended, Standard 17.3 cm


Gracey curettes are area specific curettes for the removal of deep subgingival calculus and for the finishing of subgingival root surfaces. The blade has two sharp edges.

Item	Size
1 2801071	#1-2 Incisors - All surfaces
2 2801073	#3-4 Incisors - All surfaces
3 2801074	#7-8 Premolars & Molars - Buccal lingual side
4 2801070	#11-12 Premolars & Molars - Mesial side
5 2801072	#13-14 Premolars & Molars - Distal side



Elevator Bein

Used to loosen the tooth from the periodontal ligament and ease extraction. Due to the thicker handle, this model fits particularly well in the hand.



Item	Description	Size
1 2801084	Straight 3 mm	13.5 cm
2 2801085	Round 2 mm	14 cm
3 2801086	Round 3 mm	14 cm
4 2801087	Round 4 mm	14 cm

Elevator Cryer

Designed to loosen and elevate the tooth from the periodontal ligament to make extractions easy. "Cryer-White" H-Type C. Made of Stainless steel.



Item	Description	Size
2801089	Elevator Cryer Right	13 cm

Elevator Heidbrink

Elevator Heidbrink is a delicate instrument for the removal of root tips in difficult-to-reach areas.



Item	Description
2801092	Elevator Heidbrink

Elevator Flohr

Elevator Flohr has a curved end and is used for loosening the tooth roots from the periodontal ligament and ease extraction. Available in curved-left or curved-right.



Item	Description	Size
2801091	Elevator Flohr Curved Right	15.5 cm

Root Elevator Hylin

Elevator used for the extraction of a tooth or retained roots.



Item	Description
2801247	Root Elevator Hylin 3 mm

Gingivectomy Scissor Goldmann-Fox

Scissors with one serrated (toothed) blade used for trimming tissue or cutting sutures. Used in dental procedures to trim the gingival tissues.



Item	Description
2801266	Curved Toothed 13 cm

Periosteal elevators are mainly used to lift full thickness soft tissue flaps. The tips are very sharp.

Periosteal Dental Elevator #A2



Item	Description	Size
2801234	#A2 Comfort	17.1 cm

Periosteal Elevator Moults Double-ended #9



Item	Description	Size
2801235	Moults Double-ended #9 Comfort	17.1 cm

Periosteal Elevator Double-ended Moults Octagonal Handle



Item	Description	Size
2801237	Double-ended Moults Octagonal Handle	18 cm

Periosteal Elevator Raspatory Williger



Item	Description	Size
2801236	Raspatory Williger 3 mm	16 cm

Periosteal Elevator Raspatory Freer

Mainly used to lift full thickness soft tissue flaps.



Item	Description	Size
2801241	Raspatory Freer Sharp/blunt Curved 6 / 5 mm	20 cm

Luxator



Item	Description	Size
2801177	Remeus/Fahrenkrug for Rodent Molars	18 cm

Extraction Tooth Fahrenkrug

The Fahrenkrug tooth extractor has two oppositely bent and ground lifter blades with a concave or convex cutting direction.




Item	Description	Size
2801347	For Cats and Small Dogs	19 cm

Extracting Forceps – Incisors & Canines



Item	Description	Size
1 2801101	1100 #1	17 cm
2 2801102	1100 #2	16.5 cm
3 2801108	Small #37	14.5 cm

Extracting Forceps – Premolars



Item	Description	Size
1 2801103	1100 #7 Premolars	17 cm
2 2801105	500 #13 Premolars	15.5 cm

Extracting Forceps – Molars



Item	Description	Size
1 2801106	Angled 500 #22 Molars	15.5 cm

Extracting Forceps – Roots

Item	Description	Size
1 2801104	400 #51A Roots	17.5 cm
2 2801107	500 #51A Roots	18.5 cm
3 2801112	Large Dogs	17 cm
4 2801113	Small Dogs	12 cm



Splinter Forceps Stieglitz

For removing obstructions or debris from the root canal, and removal of root fragments. Available with straight or angled tips (45 degrees).


Item	Description	Size
1 2801310	Curved	13.5 cm
2 2801311	Straight	13.5 cm



Excavator Double-ended Octagonal Handle Solid Grooved

Excavators are often used to remove small root residues and other dead tissue. Excavators are smaller than the extraction forceps and elevators and therefore easier in some situations. Excavators have a handle and often a curved shaft on both sides with a small sharp circular worktop.

Item	Description	Size
1 2801098	#127/128 1.7 mm	17.4 cm
2 2801099	#129/130 1.5 mm	17.4 cm
3 2801100	#131/132 1.2 mm	17.4 cm



Molar Cutter



Item	Size
2801178	18 cm

Bone Rongeur Forceps



Item	Description	Size
2801246	Luer Slim Curved	15 cm

Bone Rongeur Forceps



Item	Description	Size
2801065	Mini Friedman	12.5 cm

Parodontology

Hemostatic Solution

Haemostatic solution is a topical astringent solution of Aluminium chloride. Gingival retraction can be achieved using a plain gingival retraction cord moistened with Haemostatic solution. For use to aid in cessation of gingival bleeding.



Item	Description	Size
2800811	25% Aluminium Chloride	30 ml

Plaque Indicators

Plaque is a major cause of periodontal conditions but remains invisible until stained. Plaque indicators color affected areas red, making them clearly visible. This helps inform and motivate pet owners and ensures thorough plaque removal during scaling procedures. The product is simple and hygienic to use.



Item	Description	Quantity
2800862	Swabs	25 pcs

Disposable Propy Angles Universal

Our Propy Angles provide a gear design which delivers a continuous back and forth oscillation of the non-latex propy cup. In a veterinary office the benefits obtained from oscillating motion over centrifugal motion, include no costly motor repairs arising from entrapped hair and no wasted time taking measures to prevent hair entanglement. Each disposable Propy Angle comes individually wrapped.



Item	Quantity
2800791	100 pcs

Propy Paste Cups

Fluoride-free paste with fine, medium, or coarse grit for cleaning and stain removal. Non-splatter formula with unique polishing agents. For professional use only.



Item	Description	Quantity
2800819	Medium Grit	200 pcs

Hedström Files

Hedström files are being used to widen occlusal third of canal and to smooth walls of preparation.

Hedström Files Stainless Steel – 60 mm (VET)

These stainless steel Hedström files are specifically made for Animal Health and have a length of 60 mm. Assorted packs with 3 different sizes, one each.



Item	ISO	Quantity
2800730	30 35 40	Set 3 / pk
2800731	45 50 55	Set 3 / pk
2800732	60 70 80	Set 3 / pk
2800733	90 100 110	Set 3 / pk

Barbed Broaches

Set consisting of 8 nerve extrusion needles.

Use: insert the nerve needle, rotate 180° and then remove it. Working length 47 mm + handles.
Set with 4 sizes (#2-3-4-5), 2 x each.



Item	ISO	Quantity
2800734	2 3 4 5	Set 8 / pk

Application & Modelling Tools

Cement Spatula



Item	Description	Size
2801067	Octagonal Handle Solid Serrated	6.5 mm

Ball Burnisher Amalgam

A burnisher is a surgical instrument for the purpose of smoothing and/or polishing amalgam fillings and polishing composite fillings, as well as emphasizing grooves or removing excess filling material on the surface of the tooth.



Item	Description	Size
2801053	Double-ended #25/31	17 cm

Tamponade tamper



Item	Description
2801324	Luniatschek #13,0mm Comfort

Spatula Heidemann



Item	Description
2801307	#12.0/2.0 mm Octagonal Handle Solid Grooved

Miscellaneous

Mouth Spatula Steel



Item	Size
2801187	16/22 mm 17 cm

Tongue Depressors

Non-sterile wooden tongue depressors made from flexible and break-proof birch wood with smooth surface.



Item	Description	Quantity
2800167	Non-Sterile 15 cm	500 / box



Pro Pet Solutions

Inspired by Veterinarians

Pro Pet Solutions enables veterinary professionals to offer pet owners affordable, veterinarian-inspired home care products.

Veterinarians and pet owners can trust that pets receive appropriate at-home treatment. Pro Pet Solutions supports clinic profitability by encouraging repeat purchases and building customer loyalty, while ensuring long-term patient care. Our range includes flea combs, protective boots, and many more. Check out by yourself!

“Covetrus Pro Pet Solutions is the newest addition to our portfolio providing veterinarians with retail products for professional at-home care through pet owners.”

Pro Pet Solutions Index

Covetrus Pro Pet Solutions		109
1	Nutraceuticals	110
	Joint support	110
	Urinary support	111
	Skin & Coat support	111
	Immune support	112
	Liver support	113
	Stress support	113
	Puppy support	114
	Hairball Remedy	114
	Kidney support	114
	Intestinal support	115
	Rehydration support	115
	Gastro support	116
	Metabolism support	116
Dental support	116	
Other dental chew	116	
2	Topicals and Shampoo bars	117
	Topicals	117
	Shampoo bars	118
3	Patient Care & Animal Handling	119
	Protective boots	119
	Cat Litter	119
	Muzzles	120
	Ectoparasites	120



Nutraceuticals

Our NutriCareVet line comprises over 20 high-quality nutritional supplements for pets and farm animals, developed in collaboration with veterinarians to support animal health and welfare across a wide range of situations. Manufactured in Europe to the highest quality standards and based on the latest scientific research, NutriCareVet products deliver reliable quality at a fair price for pet owners and farmers alike.



NutriCareVet - Joint support

Canine Supplement - Joint Support Chewable Tablets

Joint Support chewable tablets promote musculoskeletal health and activity in dogs, especially seniors.

- It contains glucosamine, chondroitin sulfate, L-methionine, devil's claw extract, green-lipped mussel, vitamin C, and vitamin E
- These ingredients support joint metabolism, cartilage regeneration, and connective tissue development
- The formula includes 100% poultry liver, omega-3 fatty acids (EPA + DHA)



Item	Content	Net Weight
2802083	85 Chewable tablets	213 g

Canine Supplement - Joint Support Soft Chews

Soft chews for dogs and puppies (from 12 weeks) to support joint health, cartilage, and connective tissue. A concentrated blend of active ingredients that maintain healthy bones and flexible joints and helps keep active pets mobile.

- The formulation includes glucosamine, chondroitin, and egg shell membrane
- Highly palatable chicken flavour, palm oil-free, lactose-free, and free from dyes and flavour enhancers



Item	Content	Net Weight
2802245	60 Soft chews	300 g

Feline & Canine Supplement - Joint Support Liquid

Joint Support Liquid is for dogs and cats and promotes long-term joint health. Recommended for pets with stiffness, limited mobility, or lameness, especially older animals.

- It contains glucosamine, MSM, chondroitin sulfate, fish oil, vitamin C, and vitamin E
- The formula includes omega-3 fatty acids (EPA + DHA), is caramel-flavored, lactose- and gluten-free, clinically tested



Item	Content	Net Weight
2802085	Liquid	250 ml



NutriCareVet - Urinary support

Canine Supplement - Urinary Support Soft Chews

Soft chews for dogs and puppies (from 12 weeks) to help maintain healthy urinary tract, kidney and bladder function and support metabolism and digestion.

- Features Quad Action with dandelion root, nettle leaf, cranberry, and D-mannose.



Item	Content	Net Weight
2802248	60 Soft chews	300 g

Feline & Canine Supplement - Urinary Support Paste

Urinary support paste to support the physiological function of the urinary tract in dogs and cats.

- It contains methionine, taurine, vitamin E, and cranberry extract.



Item	Content	Net Weight
2802100	Paste	100 g



NutriCareVet - Skin & Coat support

Canine Supplement - Skin and Coat Support Soft Chews

Soft chews for dogs and puppies (from 12 weeks) are developed to support dry and itchy skin.

- Triple Action of omega oils – Salmon oil, Krill oil and Seaweed meal
- The chews provide essential nutrients that nourish the skin and coat which support soft and supple skin, healthy glossy coat, skin barrier integrity and replenish dry skin
- Highly palatable chicken flavor, palm oil-free, lactose-free, and free from dyes and flavor enhancers



Item	Content	Net Weight
2802246	60 Soft chews	300 g

Feline & Canine Supplement - Skin and Coat Support Liquid

Skin & Coat Support Liquid is for dogs and cats and promotes healthy skin and a shiny coat.

- It helps correct fatty acid imbalances in the diet by providing essential omega-3 fatty acids which support skin, coat, nails, and hair
- The formula includes fish oil, evening primrose oil, anise oil, and vitamins A, D3, and E



Item	Content	Net Weight
2802080	Liquid	250 ml

NutriCareVet - Immune support

Feline & Canine Supplement - Pet Multivitamins

Pet Multivitamins contain an essential blend of vitamins, minerals and trace elements for both cats and dogs.

- They are ideal for cats and dogs in stressful situations with high nutritional requirements, such as in recovery phases or when nursing
- Additionally, the tablets help maintain a healthy and effective immune system



Item	Content	Net Weight
2800629	Tablets	50 x 1g

Feline Supplement - Immune Support Lysine Gel

NutriCareVet Feline Immune Support Lysine Gel helps strengthen the immune system of cats during times of stress.

- It contains L-lysine, an essential amino acid that cats cannot produce themselves and is especially recommended for cats with herpes virus infections
- The lysine gel is palatable with no bitter taste, lactose- and gluten-free, hypoallergenic, and made from organically sourced ingredients



Item	Content	Net Weight
2802081	Gel	50 ml

Feline Supplement - Immune Support Chewable Tablets

These tablets are developed for cats of all life stages to support a healthy immune system.

- The formula contains beta-glucans, colostrum powder, L-carnitine, prebiotics, probiotics (Enterococcus faecium), taurine, and essential vitamins and trace elements
- It supports immunity during recovery, vaccinations, oral health, and weaning



Item	Content	Net Weight
2802101	190 Chewable tablets	142,5 g

Canine Supplement - Immune Support Chewable Tablets

Chewable tablets for dogs of all life stages to support a healthy immune system.

- It contains beta-glucans, prebiotics, vitamin E, vitamin C, methionine, carnitine, and Echinacea purpurea (coneflower)
- It supports immunity during recovery, vaccinations, weaning, and oral health



Item	Content	Net Weight
2802102	85 Chewable tablets	213 g

NutriCareVet - Liver support

Feline Supplement - Liver Support Chewable Tablets

NutriCareVet Feline Liver Support chewable tablets promote healthy liver metabolism in cats.

- It contains milk thistle, artichoke, betaine hydrochloride, vitamins E and B12, arginine, and taurine
- These ingredients support liver function, digestion, and overall vitality



Item	Content	Net Weight
2802087	80 Chewable tablets	60 g

Canine Supplement - Liver Support Chewable Tablets

NutriCareVet Canine Liver Support chewable tablets promote healthy liver metabolism in dogs.

- It contains milk thistle, turmeric, artichoke, betaine hydrochloride, lecithin, schisandra, and vitamin E.
- These ingredients offer antioxidant, anti-inflammatory, and regenerative benefits, support bile function, stimulate metabolism, aid digestion, and protect the liver



Item	Content	Net Weight
2802088	85 Chewable tablets	213 g

NutriCareVet - Stress support

Feline and Canine Supplement - Stress Support Liquid

NutriCareVet Stress Support Liquid is developed for dogs and cats to help reduce anxiety in stressful situations like thunderstorms, fireworks, loud noises, or travel.

- It contains L-tryptophan and green tea extract to support serotonin production, promoting calm and relaxation, along with vitamin B1 to help support nervous disorders
- The caramel-flavored formula is non-sedating, hypoallergenic, lactose- and gluten-free, dye-free, and clinically tested
- It is easy to administer and effective within 90 minutes



Item	Content	Net Weight
2802074	Liquid	250 ml

Feline and Canine Supplement - Stress Support Tablets

NutriCareVet Stress Support Tablets are developed for dogs and cats that help calm anxious or nervous pets during stressful situations such as thunderstorms, fireworks, loud noises, or travel.

- The tablets contain tryptophan, taurine, and vitamins B3, B6, and D3 to support emotional balance and relaxation
- They are non-sedating, lactose- and gluten-free, dye-free, clinically tested



Item	Content	Net Weight
2802073	Tablets	43 g

NutriCareVet - Puppy support

Milk Replacer Powder

A complete supplement for young or weakened puppies and convalescent adult dogs.

- Supports intestinal health and immune function
- Closely resembles the composition of natural breast milk
- It contains whey protein, vegetable oils, probiotics, prebiotics, vitamins and trace elements



Ideal bundle opportunity!

Nursing Kit for kitten and puppies

Ergonomically shaped mouthpiece. The top of the bottle can be easily removed thanks to the screw cap and is therefore ideal for cleaning with the brush supplied.



Item	Content	Net Weight
2802625	Powder	340 g

Item	Quantity
2802563	30 pcs

NutriCareVet - Hairball Remedy

Feline Supplement - CatLax Hairball Remover

CatLax paste helps the flow of ingested fur through the gut to prevent hairball formation.

Cats can easily pick up hair from frequent grooming, which causes hairballs to form. Ingested hair can become lodged in the digestive tract and form a cluster or hairball. Symptoms of hairballs include choking, coughing, vomiting, and constipation. Hairballs can be avoided with frequent grooming and routine use of CatLax.

- CatLax contains ingredients that gently lubricate the digestive tract to prevent hairball build-up
- The paste acts as a lubricant that coats the swallowed hair and allows the hairball to travel down the digestive tract
- It has a mild laxative effect that penetrates hardened stool so it can be passed



Item	Content	Net Weight
2802450	Paste	60 g

NutriCareVet - Kidney support

Feline Supplement - Charcoal Paste

Supplement for aging cats that helps bind harmful metabolites in the lower gastrointestinal tract and support kidney function.

- The paste contains activated charcoal (50%) and glycerin
- It helps bind uremic toxins, which are produced as byproducts of bacterial fermentation in cats' intestines
- These toxins are absorbed by the charcoal and excreted in the feces, rather than being absorbed by the intestines and excreted via the kidneys



Item	Content	Net Weight
2802665	Paste	30 ml

NutriCareVet - Intestinal support

Canine & Feline Supplement - Pet Intestinal Formula Paste

Pet Intestinal Formula is a complementary mineral paste developed for dogs and cats that supports intestinal flora after digestive disorders.

- The immediate response to restore normal digestive function
- The paste helps to decrease volume and frequency of stool output
- It restores electrolyte and water balance while provides a source of easily assimilated energy
- It helps repair and restore gastrointestinal mucosa
- It contains bentonite-montmorillonite, prebiotics, probiotics, electrolytes

Item	Content	Net Weight
2800946	Paste	15 ml
2800947	Paste	30 ml
2800949	Paste	60 ml



Canine Supplement - Intestinal Support Soft Chews

Soft chews for dogs and puppies (from 12 weeks) to support a healthy gastrointestinal system, gut flora, and digestion.

- It contains psyllium husk, prebiotics, probiotics, B vitamins, vitamin E, salts, and minerals
- It includes an enzyme blend to aid digestion of starches and cellulose
- Over 1 billion probiotic Enterococcus faecium bacteria per chew
- Fast-acting support for dogs with sensitive stomachs

Item	Content	Net Weight
2802247	60 Soft chews	300 g



NutriCareVet - Rehydration support

Canine & Feline Supplement - Oral Rehydration Powder Sachets

Oral Rehydration Powder Sachets for dogs and cats is designed to restore digestive health and support recovery.

- It provides essential vitamins, minerals, and amino acids, making it ideal for animals with digestive disorders or in need of intensive care
- Free from citrate, it supports the reduction of acute intestinal issues and is suitable for pets of all life stages, breeds, and sizes
- Simple to administer and enriched with B vitamins to normalize metabolism
- It is a complete and balanced solution for rehydration and recovery

Item	Content	Net Weight
2800950	Powder	10 x 30 g



Bovine Supplement - Rehydration Support Tabs

This oral solution is developed for calves to treat dehydration caused by scours.

- It contains prebiotics, probiotics, sodium bicarbonate, and lactose to stabilize water and electrolyte balance and restore energy
- Easily soluble tablets with a SID value of 82 mEq/litre
- Individually wrapped for simple application

Item	Content	Net Weight
2802096	30 Tabs	1500 g



NutriCareVet - Gastro support

Feline Supplement - Gastro Support Chewable Tablets

Chewable tablets for cats to support gastric function and digestion, especially in cases of acidosis.

- With psyllium husks, elm bark, marshmallow and montmorillonite
- Stimulation of the bowel movement through probiotics and prebiotics

Item	Content	Net Weight
2802092	80 Chewable tablets	60 g



Canine Supplement - Gastro Support Chewable Tablets

Chewable tablets for dogs to support gastric function and digestion, especially in cases of acidosis.

- With psyllium husks, elm bark, marshmallow and montmorillonite
- Stimulation of the bowel movement through probiotics and prebiotics

Item	Content	Net Weight
2802093	85 Chewable tablets	213 g



NutriCareVet - Metabolism support

Bovine Supplement - Metabolism Support Bolus

Complementary feed for cows to support metabolism and calcium levels post-calving, when calcium demand is increased. Contains calcium formate and vitamin D3.

Item	Content	Net Weight
2802105	8 Boli	approx. 880 g



NutriCareVet - Dental support

Canine Supplement - Dental Support Soft Chews

Soft chews for dogs and puppies (from 12 weeks) to support healthy teeth, gums, and fresh breath.

- It contains omega-3 and -6 fatty acids, seaweed meal, spirulina, cranberry, pomegranate, green tea extract, and various oils
- The chewy texture and abrasive ingredients help reduce plaque
- Highly palatable chicken flavor, palm oil-free, lactose-free, and free from dyes and flavor enhancers

Item	Content	Net Weight
2802249	60 Soft chews	300 g



Other Dental Chew

Ostrich Twister Dental Chew

100% all-natural dog chews made from the tendons of South African Ostriches. Ostrich Twister chews do not gelatinise like other chews making them extremely palatable. As the dog chews, the fibres from the tendon move around the teeth creating a natural flossing effect which makes Ostrich Twisters the ideal dental chew.



Item	Size
2800745	10 pcs

Topicals and Shampoo bars

Topicals

Ear Cleaner with Flexi Nozzle

Gentle, non-toxic ear cleaning solution for dogs and cats. Removes dirt and wax, keeps ears clean and dry, and helps prevent staining around the ears.

- Easy to use with the special flexible nozzle
- Limits any staining or discolouring of hair around the ears
- With Aloe Vera



Item	Net Weight
2800630	100 ml

Single-Use Eye Cleaning Wipes

Eye cleaning wipes are very soft and are impregnated with a fragrance-free solution which moisturizes the eyelids and eyelashes.

- The wipes can be used
- at room temperature,
 - refrigerated (for instance for swelling) or
 - can be heated in its packaging in the microwave (up to 40 °C) for an even more effective cleaning.



Item	Content	Quantity
2802433	Wipes	20 pcs

Eye Area Cleaner Hydrogel

For use around the eye area to clean and soothe.

- It effectively removes tear stains and debris
- Is alcohol-free, pH neutral, non-irritating, and a safe, effective solution



Item	Net Weight
2800951	60 ml

Eye Drops

Preservative-free drops with soothing, refreshing properties for dry, red, or irritated eyes. Suitable for all animals.

- Made with 100% natural ingredients including chamomile, euphrasia, and hamamelis



Item	Net Weight
2802431	10 ml

Paw and Nose Balm

Paw and nose balm forms a protective barrier on pets' pads and increases their resistance to friction. Paws can be dry and damaged by unprotected contact with elements such as ice, salt, sand, rough, hot surfaces or after long walks.

- The balm is organic, vegan, 100% natural vegetable blend of butters, special plant oils and dragon's blood
- Without petroleum derivatives
- Non-toxic and harmless when licking
- It nourishes and moisturizes rough, cracked and dry paws and nose gently



Item	Content	Net Weight
2802432	Balm	50 ml

Shampoo bars

Eco friendly Pet Shampoo Bars

- Made from natural and naturally-derived ingredients
- pH-balanced for a clean and gentle wash
- Easy to use with less rinsing needed
- No plastic packaging
- 100% vegan
- Eco-friendly
- Degradable without residue
- Suitable for all pets



Pet Shampoo Bar - Sensitive scent-free

- Ideal for sensitive skin and noses
- Contains natural coconut oil and panthenol



Item	Description
2802426	Sensitive scent-free

Pet Shampoo Bar - Keep-insects-away with lemon scent

- Ideal to help repel ticks and fleas
- Contains citronella for its naturally bug-repellant properties and fresh scent, as well as Neem which is anti-fungal and kind to the skin



Item	Description
2802427	Keep-insects-away with lemon scent

Pet Shampoo Bar - Relaxing with lavender scent

- Ideal to help relax while washing
- Contains calming lavender and hibiscus petals to soften and add shine to the coat and make wash times more relaxing for your pet



Item	Description
2802428	Relaxing with lavender scent

Pet Shampoo Bar - Odour-neutralising with mint scent

- Ideal to neutralize odours
- Contains tomato extract and mint to neutralize odors



Item	Description
2802429	Odour-neutralising with mint scent



Patient Care & Animal Handling

Protective boots

Petcare Paw

Item	Size
2800821	130 x 60 mm, S
2800822	150 x 75 mm, M
2800823	177 x 85 mm, L
2800824	215 x 90 mm, XL
2800825	230 x 100 mm, XXL
2800826	230 x 110 mm, XXXL



Size	Breed
S	Yorkshire Terrier, Toy Poodle
M	Jack Russel Terrier, West Highland Terrier
L	Border Collie, Beagle, Dalmatian
XL	Labrador, Boxer, Setter
XXL	Rottweiler, Great Dane, German Shepherd
XXXL	Newfoundlander, St. Bernard, Leonberger

Cat Litter

Cat Litter Bentonite

The cat litter is a natural clay (bentonite). Forms solid and compact lumps when wet.



Item	Description	Size
2800631	Clumping Light-grey	8 kg

Silica Crystals

- Biodegradable
- Antibacterial
- Composition: 100% silicon dioxide
- Ultra-absorbent



Item	Description	Size
2800883	Lumpy Silica	2 x 16 L
2800884	Lumpy Silica	5 L
2800885	Micro Silica	5 L

Recycled Paper Cat Litter

- Made from 100% recycled ink free paper
- Biodegradable
- High odour control
- Exceptionally absorbent

Item	Description	Size
2803129	Recycled Paper Cat Litter	30 L



Muzzles

Dog Muzzle for Short-Nosed Breeds

- Adjustable snap-lock clip closure
- Easy to fit securely



Item	Width (cm)	Height (cm)	Size
2802598	20.3	9.8	S
2802599	22.5	11	M
2802600	24.6	13.5	L
2802601	29	15	XL

Dog Muzzle Nylon

- Adjustable snap-lock clip closure
- Easy to fit securely

Item	Size
2802305	N.0
2802306	N.1
2802307	N.2
2802308	N.3
2802309	N.4
2802310	N.5
2802311	N.6/4XL
2802312	N.7/5XL
2802313	N.8 Short



Ectoparasites

Fleacomb



Item	Description
2800665	Without Handle
2800666	With Handle

Tick Remover Pen-Style Plastic



Item	Description
2800667	Pen-Style Tick Remover

Tick Forceps Crossed Steel



Item	Size
2801325	14.5 cm

TICKLESS - Ultrasonic Tick and Flea Repellent

Chemical and fragrance free tick repellent. Works with unaudible soundwaves which irritate ticks and cause them to fall off the pet. Provides at least 6 months of protection and is safe for animals of every age.

Item	Quantity
2802749	1 pc



NutriCareVet

from Covetrus

SUPPLEMENTS AS GOOD FOR YOUR PETS AS THEY TASTE!

- High quality ingredients & affordable price
- Highly palatable and suitable for puppies and senior pets
- Developed with and recommended by veterinarians
- Made in Europe by leading manufacturers
- Sustainable packaging
- Easy support for optimal health from home
- Made in Europe by leading manufacturers



OUR ASSORTMENT

Joint	XX		X	
Immune	X	XX	X	
Intestinal	X		XXX	
Rehydration			X	X
Gastro	X	X		
Skin and coat	X		X	
Kidney		X		
Liver	X	X		
Hairball Remedy		X		
Dental	X			
Stress			XX	
Urinary	X		X	
Metabolism				X
Puppy	X			

CARE WITHOUT COMPROMISE

